

Product datasheet for **MC227166**

Fam71f1 (NM_001289663) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam71f1 (NM_001289663) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam71f1
Synonyms: GARI-L1; Gm764; NYD-SP18
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227166 representing NM_001289663
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGATGACATCAGTTCACCTAGAAAGTCTTGGTGGACATCGAAGAAGACAGTAAAAGTGATAAGATCAT
ATCCAACCTTCCTTCCCTTGAATGCCTGGGAGAAGTTCAGGGGCTCTTGCTGTGGATGGGAAACGAA
CCCTGGAGTGGGCTGGGAGTGGAGGAAGGACTACTCTGCCAGATGGTTCATTCTCCAGAATCAATCTG
TTTCCTAACTCCGTGGTGTTCGAAAGCAACTTTGTTCCAGTCCAGAGGAGCAGAACTGGAAAGAAATCT
ATAAGGCTTCCACACAATGGCCCTTGGGGTGACCTCCTCCGTACCTTGCTTCCCTTCCGAACATCCT
CCTTATGGCCAGAGTCAAATGGCATCAGGGGCAGAGCCAGACATGGAACAGACCATCCAGAGCGCCAAGC
ATCAACCTGAAGAGCATTCTCCCTTTGAAGTTTGTGGAGCTCCAGATCTGGGACGACCAGGAACGCGTCC
TGCGGTTGAGGACAGTCAACGAGAAAATCTACTATCTGAAGCTCCATCCCGACCACCCTGAGACTGTCTT
CCACTTTTGGATCCGCTGGTCAAATCTGCATAAGGGTCTGTCCATCACCACCAAGGACCCTACCATC
CTTGTCACCTCACTGCCTGGTACCAAGAGCCTTTGCAGCCAGGTGAAAGACTGAGTTAGTGCAGAAGA
AATCGAAAGGCCTCCAGCCAGTGAGAGCCTCACACATCTGATGGCCCAGGGGAGAGCCAGCCCTCTC
CCAGATCTTTTCCGACTTGCATCAGCAGAAAACAGTACAGTTCCAGCAGGAGCGAAAAGATGCACATCAAC
AAGACCAGCTCAGAGAAAGCTACTCCCTGTGAAGACAGCATCCCCTGTACCTGTGATCTCAACTGGAGGG
ATGCATTATGTTTCGGAGAATGGGAAAGAGAAAACCTTCTGGACCACAGCCCCTCTCACTCCTCGGTAC
TCTGGCTGCCTCCAGCAGGCCACGGCTGTCCCTAACTGGAGGAACTCTAT**TGA**

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_001289663



[View online »](#)

Insert Size:	1035 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289663.1</u> , <u>NP_001276592.1</u>
RefSeq Size:	1596 bp
RefSeq ORF:	1035 bp
Locus ID:	330277
UniProt ID:	<u>Q3UZD7</u>
Cytogenetics:	6 A3.3