

Product datasheet for MC227220

Rcor3 (NM_001290279) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rcor3 (NM_001290279) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rcor3
Synonyms: 4921514E24Rik; C730034D20Rik; E130101E15Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227220 representing NM_001290279
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCCGCATGATGGAGAAAGGGCCCGAGTTACTGGGAAGAGCCGATCGGCCAACGGCGGCCAAGA
 GCCCGGGGGAGGCGGGCTCCTCGGCCAACGGCGGGCTGCACTTCTCGGAGCCGAGAGCGGCTGCAG
 CAGCGACGACGACGACGGGGATGTTGGGATGAGAGTCGGAGCTGAATACCAAGCTCGGATCCCTGAATTC
 GATCCAGGTGCTACCAAGTACACAGATAAAGACAATGGAGGCATGCTTGTGTGGTCTCCTTATCACAGTA
 TCCAGATGCCAAATTGGATGAATATATCGCCATTGCAAAAGAGAAACATGGCTACAATGTTGAACAGGC
 GCTTGGCATGTTGTTCTGGCATAAACACAACATTGAGAAGTCCCTTGGCGACCTCCCTAACTTACGCCCC
 TTTCCGGATGAGTGGACAGTGGAAAGACAAGTCCCTGTTTGAACAAGCCTTTAGTTTTTATGGGAAGAGCT
 TTCACAGGATTCAGCAAATGCTTCCAGATAAGACAATTGCAAGTCTTGTAAAGTATTACTATTCTTGAA
 AAAAATCGCTCAAGGACAAGTTTGATGGATCGCCAGGCTCGTAACTAGCCAATAGACACAATCAGGGT
 GACAGTGATGATGACGTAGAAGAAGCTCACCCATGGATGGAAACGACAGTGACTACGACCCAAAGAAAG
 AAGCCAAAAGAGAGGGTAACGCTGACCAGCCGGTCCAAACCAGCAAGATCGGCCTGGGCGGAGAGAGTA
 CCAGAGTCTGCAGCACCGCCACCACTCTCAGCGCTCCAAGTGCCGCCCTCCTAAGGGCATGTATCTAACCC
 CAGGAGGACGTGGTGGCTGTCTCCTGCAGCCCCAACGCAGCCAACACCATCCTGAGGCAACTGGACATGG
 AGCTGATCTCTCAAGCGCCAGGTTGAGAATGCTAAGCAAGTAAACAGTGCCTTAAGCAGAAAATGGA
 AGGTGGAATTGAAGAGTTCAAACCGCCTGAGGTGTCGCAAAATATGGTAAAGATTTTCAAGCTATTGCAG
 ATG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001290279



[View online »](#)

Insert Size:	1056 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001290279.1, NP_001277208.1</u>
RefSeq Size:	1637 bp
RefSeq ORF:	1056 bp
Locus ID:	214742
Cytogenetics:	1 H6
Gene Summary:	May act as a component of a corepressor complex that represses transcription. [UniProtKB/Swiss-Prot Function]