

Product datasheet for **MC227235**

Gas7 (NM_001277080) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gas7 (NM_001277080) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gas7
Synonyms: AW124766; B230343A10Rik; Gas7-cb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227235 representing NM_001277080
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCAACACGACTACCAATGAGACCCTGGGAACGTCCAGCAGTTCTCCTGGGATTCAGCCAGCCC
TGGCCCTCACAGGAGCTCTTTGCCTACAACAGAACCTGGGATCCTCATCGCCAGGCAGGAAACAGAGCAA
GGAAAACACCATCACCATCAACTGTGTGACATTCCTCACCCGGACACGATGCCTGAGCAACAGCTGCTC
AAGCCACCGAGTGGAGCTACTGTGACTACTTTTGGGCGGACAAGAAGGCCACAAAGGCAATGGCACGG
TGCTGGCTTTGAACTGCTGCTGCAAGAACAAGTAAAGGCAAGCAGATGCAGAAGGAGATGTCGGAGTT
CATCCGGGAAAGGATAAAGATTGAAGAAGAATATGCAAAGAACCTGGCTAAGCTCTCTCAGAACTCCTTG
GCTGCCAGGAGGAAGGCTCCTTGGGAGAGGCATGGGCCAAGTAAAGAAGAGCCTGGCAGATGAGGCTG
AAGTTCATCTCAAGTCTCTGCCAAGCTCCACAGCGAGGTGGAGAAACCCCTGATGAACTTCAGAGAGAA
CTTCAAGAAAGACATGAAAAAGTGCACCACCACATCGCTGACCTCCGTAAGCAGCTGGCAAGCCGCTAT
GCCTCAGTGGAGAAGGCCCGCAAAGCCCTCACAGAGCGGCAGAAAGACCTGGAGATGAAGACCCAGCAGC
TGGAGATCAAGTTGAGCAACAAGACCCGAGGAGGACATCAAGAAGGCACGGAGGAAGTCCACTCAGGCCGG
AGATGACCTCATGCGCTGTGTGGACCTCTACAACCAGGCCAGTCCAAGTGGTTTGAAGAGATGGTGACG
ACCACACTGGAAGTGGAGCGGCTGGAGGTGGAGAGGGTGGAAATGATCCGACAACATCTGTGCCAGTACA
CACAGCTGCGGCACGAGACGGACATGTTCAACCAAGCACAGTCGAGCCTGTGGACCAACTGCTGCGAAA
AGTGGACCCAGCCAAAGACAGAGAGCTATGGGTCAGAGAGCACAAAACAGGCAACATCCGCCAGTGGAC
ATGGAGAT**TAG**

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001277080



Insert Size:	1062 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277080.1, NP_001264009.1</u>
RefSeq Size:	6848 bp
RefSeq ORF:	1062 bp
Locus ID:	14457
Cytogenetics:	11 B3