

## Product datasheet for **MC227263**

### Nipa2 (NM\_001256132) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nipa2 (NM\_001256132) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nipa2  
**Synonyms:** 2600017P10Rik; 3830408P04Rik; AB041581  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227263 representing NM\_001256132  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAGCCTGGGGCGTGGAAAATGATTTTTATATTGGTCTGGGATTGGCTATGACCTCCAGCATTTTCA  
TTGGAGGGAGTTTCATTTTGAAGGAAAGGCCTTCTGCGACTTGCAGGAAAGGCTCAATGAGAGCAGG  
TCAAGGAGGCCATGCATATCTTAAGGAATGGTTGTGGTGGCTGGACTGCTGCAATGGGAGCCGGTGAA  
GTGGCCAATTTTGTGCTGCATATGCATTTGCTCCGGCCACATTAGTGACTCCACTAGGTGCGCTCAGCGTCC  
TCGTAAGTGCCATTCTGTCTTCACTTCTAAATGAAAGACTTAATCTTCATGGGAAAATTGGGTGTT  
GCTAAGTATTCTAGGGTCTACAGTTATGGTCATTCATGCTCCAAAGGAAGAGGAGATTGAGACCTTAAAT  
GAAATGTCACATAAGTTAGGTGATCCAGGTTTTGGTCTTTGCAACATTTGTGGTCATTGTGGCCTTGA  
TATTCATCTTTGGTGGGACCTCGACATGGACAGACAAACATCTTGTGTATATAACAATCTGCTCTGT  
GATTGGAGCATTTTCAGTCTCCTGTGTGAAGGGTTGGGTATCGCCATAAAGAGTTGTTGGCTGGAAAG  
CCTGTTCTGCAGCACCCCTCGCCTGGATTCTTCTGTTTAGCCTTGTGGTCTGTGTGAGCACACAGATTA  
ATTACCTAAACCGGGCTCTGGACATTTTCAACTTCCATTGTGACTCCAATATATTATGTTCTTTAC  
AACTTCAGTTTTGACTTGTTCAGCTATTCTTTTTAAAGAATGGCAGGATATGCCTGTTGATGATGTCATT  
GGGACTTTGAGTGGCTTCTTCACAATCATTGTGGGATATTTTTGTTGCATGCCTTTAAAGATGTCAATT  
TTAGTCTAGCAAGTCTGCCAGTATCTTTTCGAAAAGATGAAAAGGCAATGAATGGCAATCTTTCTAGTAT  
GTATGAAGTTCTTAATAATAATGAAGACGACTTACCCTGTGGAATTGAACACACTGGTAAAAATATTTCC  
CGTAGGAATGAAATCTGCCATCTTTT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001256132



[View online >](#)

<b>Insert Size:</b>	1080 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001256132.1</a></u> , <u><a href="#">NP_001243061.1</a></u>
<b>RefSeq Size:</b>	3073 bp
<b>RefSeq ORF:</b>	1080 bp
<b>Locus ID:</b>	93790
<b>UniProt ID:</b>	<u><a href="#">Q9JJC8</a></u>
<b>Cytogenetics:</b>	7 B5
<b>Gene Summary:</b>	Acts as a selective Mg(2+) transporter.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Variants 1-4 encode the same isoform (1).