

Product datasheet for **MC227275**

BC024978 (NM_001243888) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: BC024978 (NM_001243888) Mouse Untagged Clone
Tag: Tag Free
Symbol: BC024978
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227275 representing NM_001243888
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACTTCTCCATGCTCTTTCCCCCTAAAACCCACAATCTCTCCAATAATTCATGAGACCCAGACACCA
ACATTCCACCTCCTTTACCTCTAAATCCCCGGACCTAGCTTTACCATCCCCTCCTGTAGTCTCCATAC
ATCTATTTTCATCCCACCTGCGCCCTCCCCCTCCTGCTCCTCCCCCTCCCCCTCCCCCTCCCCCT
CTACCCTCTGCTGTGGAGCCTGTTCTCCCTCATGTCTATGGCCTGAAGAACAGCCAGCTCCTGAAGGAGG
CCTTAGAGAAGGCTGGCCAGCCCCAAAGGCAAAGAGGATGTGAAGAGGCTCCTGAAGTGCACAAGGA
CCGGTTTCGAAGTGACCTGCAGTGGATCCTTTTTGTGCTGACCTGCCATCATGCATCCAGGAAGGCCCT
CAGTGTGGACTGGTGGCCTTGTGGATGGCAGAAGCTCTTCTGTCAACACCTGACAGCGTCTCTCTGGAAA
GACTAGTGCAGGTGGCCAAGGAGAGAGGCTACACCGCACAGGGAGAAATGTTCTCAGTGGCTGACATGGC
CAAAGTGGCCAGGAGACTTTGGACTGCCAGGCTGAGCTCCTGTGTGGTGGCTTGGTGGTCCCAACAGA
GAACGAGTTCTACAGCACCTCATCACTGGACATCCCCTGCTCATTCCCTATGATGAGGATTTCAACCAG
AGCCGTGTCAAAGAAAGGCCATAAGGCCCACTGGCAGTGAGTGCAGGGTCTGATAGGTGTGCAGAA
TGTGCCAAGCCCTGGCTACATTGAAGACTCTGAGCTACCAGGCCCTTCCACCCAGTGCCTGGTGCACCC
CACCAACCACCATCCTTCCCAGAGGAAAGCTCCCAGGTGCCCTTCTCTCTCCAAGCAGGGCAAGA
GTTGGCATTACCAGCTGTGGACTACAGCCAAGTCCGGGAGAGCAACCTGCAGCTGACAGACTTCTCACC
CGCAAGGGCTGCCGATGGACAGGTGTACGTGGTACCTGCTGGTGGGTAGAGGCTGGCCTCTGTGCCAG
GCCCTGCTGCTCAGACCAGGAGGGAAGCCAT**TAG**

ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001243888
Insert Size: 1086 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001243888.1, NP_001230817.1</u>
RefSeq Size:	1748 bp
RefSeq ORF:	1086 bp
Locus ID:	414069
Cytogenetics:	7 A3