

Product datasheet for MC227350

Gm11127 (NM_001199967) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm11127 (NM_001199967) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm11127
Synonyms: ENSMUSG00000079492
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227350 representing NM_001199967
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGGACCCTTACGAGTTGCCCTCTACTTCTCCTGCTTCTGGCGACCGCCCTGGCCCCGACCAAGCGG
 GCAAGAGCTCACACTCGCTGCGGTATTTCCGCACTGCCATATCCAGGCCAGGCTTCGGGGAGCCTCGGTT
 CACCGCGTGGGTTACGTGGACGACACAGTTTCATGCGCTTCGACAGCGACTCTGAGAATCCCAGGGCC
 GAGCCATGCAAGCCATGGGTGGAACAGATGGAGCCGAATATTGGGAGCAGGAGACGCGGAAATCAAGG
 AGCACACCCAGAATTTCCGAACCTGCCTTTACAACCTACTCCACCTGTATAACCAGAGTCAGGACGGTCC
 CCACACAATCCAGGATATGCACGGATGCTACGTGGGACCAGATGGTCAGTTCTACACGGGCACTATCAA
 CATGCCTATGATGGCCACGATTACATCACCTGAACGAGGACCTGAGCTCCTGGACCGCAGCAGATGCAG
 TGGCTCAGATCACTCAACACAAATGGGAGGAGGCTGGAGTGGCAGAGGAATACAAGGCCTACCTAGAGGG
 CACATGCGTGGAGACTCCACAGATTGTTGGAGGAAGCTCAGCAACTCAGACCCCAAAATACACAT
 GTGACCCGTACCCAGACCTGAAGTGATGTCACCTGAGGTGCTGGGCCCTGGGCTTCTACCCGTGCTG
 ACATCACCTGACCTGGCAGTTGAATGGGAGGAGCTGACCCAGGACATGGAGCTTGTGGAGCAAGGCC
 TGCAGGGGATGGAACCTTCCAGAAGTGGCAGCTGTGGTGGTGCCTCTTGGGAAGGAGCAGAATTACACA
 TGCCATGTGCACCATGAGGGGCTGCCTGAGCCACTCACCTGAGATGGGAGCCTCTTCCATCCACCGACT
 CTAACATGGTAATCATTGCTGCTCTGGTTGCTTTGGAGCTGTGATCATCATTGGAGCTGTGGTGGCTTT
 TATGATGAAGAGGAGAAACACAGGTGAAAAGGAGGGTCTATTGTTGGGAGCCGCGCCACATTCGCCG
 TTACAAGATGGCGCTGACAGCTGTGTTCTAAGTGGTAACAAATAATCTGCGCATGTGCCAAGGGTATCT
 TATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



ACCN:	NM_001199967
Insert Size:	1125 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199967.1</u> , <u>NP_001186896.1</u>
RefSeq Size:	1344 bp
RefSeq ORF:	1125 bp
Locus ID:	100529082
Cytogenetics:	17 B1