

## Product datasheet for **MC227389**

### Gm10256 (NM\_001270510) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gm10256 (NM\_001270510) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Gm10256  
**Synonyms:** ENSMUSG00000069031  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227389 representing NM\_001270510  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCAGAACTGATCAGCCTGGGAAAATTTTTATAGGAGGCCTCAACATAAAGACCAGACAAAAGACTC  
 TTCAAGAAATATTTGGGAGATTTGGACCAGTAGCACGTGTTATTTTATGAGAGATCGTGAAGAACTAAGAA  
 GTCTAGAGGCTTTGCTTTCTTACTTTCCGACGCCTTGCAGATGCTAAGAATGCTGTCAAAGAAATGAAT  
 GGAGTGATTTTGGATGGAAAAAGAAATTAAGTAAAACAAGCCAGGAGACCATCATCACTTGAAGTGGTA  
 GCAAGAAGAGACCACCATCCTTTTCAAGAACCAGAGGTGCTTCAAGAATTCTGAAATGTGGAAGAGGAGG  
 AAGAAGCAGAGCAAGAAGTGGTCCTTATGTGAAGGGAACCTGGGTGGTATAGATATACTCCTAATTTT  
 AATGTGAGTTCTTCTGGGAGACACTTTGCAGTTAAAAGAAATCCATCTTCAAAAAGGGATGATCCCCCTT  
 CTAAAAGATCTGCAACTTCTGCTCAGACAAGAAGCAACTACTGGACTAAGAGGACGAGAACCACATCGAAG  
 AGAGATTTCTAGAAATATGCCAAGAGGAGAGCCGGCGTCTTCCAGAAGAGATGAGTATCCATTACCAAGG  
 GATTATGGCCAGTCTTCTAATGATAGGAAATATGAATCAACATCCAGAGGATACTGTGACTATGGCAATT  
 ACCACTCTCGGGAAGAAAGTGCATCTAAAGTCTTCAGTGATCATGCTGGCTATCTTGGAGACGTGATAG  
 AGATTTCTCTGAATATTTAAGTGGAAATTCATACAGAGATACATACAGAAGTTATGGAAGATTCATGAG  
 GCACCATCTGCAAGAGGTGGAAACAACCCTATGATGACTATAGTAACTACAAGATGGATATGGTGGAA  
 GAGGAGAACCATACATAAGCAACCGAAGTAACATATACTCAAGTGACTATGAACGTAGTGGTAGACAAGA  
 AGTGCTTCCACCACCTATTGATAGGGAATACTTTGACCGTGAGGGCAGACAGGAAAGAGGGCATTACCT  
 AAGGATGGGCTGTATAGTGCATCTCGTGAATCATACAGTAGTAATACCAAGATTTGGGGATCCCATGGA  
 GGTCGTGGAGAAAGCAGAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_001270510
<b>Insert Size:</b>	1143 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001270510.1</a></u> , <u><a href="#">NP_001257439.1</a></u>
<b>RefSeq Size:</b>	1580 bp
<b>RefSeq ORF:</b>	1143 bp
<b>Locus ID:</b>	100042874
<b>Cytogenetics:</b>	Y