

## Product datasheet for **MC227435**

### **Nap1l4 (NM\_001285489) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nap1l4 (NM\_001285489) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nap1l4  
**Synonyms:** 2810410H14Rik; AI316776; D7Wsu30e; Nap2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227435 representing NM\_001285489  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCAGAAAACAGTCTTTCAGATGGAGGCCCTGCAGATTCTGTGGAAGCTGCTAAAAATGCCAGTAATA  
CAGAAAAACTCACAGACCAGGTGATGCAGAACCACAAGTCTGGCAGCTTTCAGGAACGTCTTGACAA  
CGTCTCTCACACTCCTTCTAGTTACATAGAACTTTACCTAAAGCTGTCAAAGAAGAATTAATGCCCTG  
AAGCAGCTCCAGGTGAGGTGTGCGCACATAGAAGCTAAGTTCTACGAGGAAGTTCATGACCTGGAGAGGA  
AGTATGCAGCCTTGTAACAGCCTCTATTTGACAAGAGAAGAGAATTCATCACTGGTGACGTGGAGCCAC  
AGATGCAGAGTCAGCGTGGCACAGTGAGAATGAGGAGGAGGACAAGTTGGCTGGAGATATGAAGAATAAA  
GTAGTCATAGCTGAGAAAGAAGCAGCAACAGTGGAAAGAGCTGAACCCAAAGGCATCCCTGAGTTCTGGT  
TCACCATCTTCAGAAATGTGGACATGCTTAGTGAGCTGGTGCAGGAATATGATGAACCGATCTTGAAACA  
TCTGCAGGATATTAAGTGAAGTTTTTCAGACCCCTGGACAACCTATGTCTTTTGTGCTAGAGTTTCACTTT  
GAACCCAAAGCAGTACTTCACCAATCCTGTCCTGACGAAAACGTACAAGATGAAATCAGAGCCAGACAAGG  
CCGACCCGTTTTCTTTGAAGTCTGAAATTGTGGACTGTGATGGGTGTACAATAGACTGGAAGAAGAG  
AAAAAATGTCACAGTCAAACCATCAAGAAAAAGCAGAAGCACAAGGGCCGGGGCACCGTACGAACCATC  
ACCAAGCAAGTGCCCAATGAGTCCTTCTCAACTTCTTCAAGCCCTGAAAGCCTCTGGAGATGGAGAAT  
CTCTGGATGAAGATCCGAGTTCACCTTAGCCTCTGACTTTGAAATTGGGCATTTTTCTGAGCGCAT  
CGTGCCACGGGTGTCTCTACTTCACTGGGGAGGCCATAGAGGATGATGACAATTTTGAAGAAGGCGAG  
GAAGGTGAAGAGGAGGAATTGGAAGGTGATGAAGAAGGAGAAGATGAGGATGATGCTGATGTTAACCCCA  
AGAAAAGAACCCAGCCAGCCAGCGGAGTGAAGCAGCAGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_001285489
<b>Insert Size:</b>	1161 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001285489.1</a></u> , <u><a href="#">NP_001272418.1</a></u>
<b>RefSeq Size:</b>	2299 bp
<b>RefSeq ORF:</b>	1161 bp
<b>Locus ID:</b>	17955
<b>UniProt ID:</b>	<u><a href="#">Q78ZA7</a></u>
<b>Cytogenetics:</b>	7 88.24 cM
<b>Gene Summary:</b>	Acts as histone chaperone in nucleosome assembly (PubMed:28366643). In condensing spermatids, mediates the loading of the heterodimer composed of histones H2AFB1 and HIST1H2BA/TH2B onto the nucleosomes, thereby promoting the replacement of histones to protamine in male germ cells (PubMed:28366643).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (a). Variants 1 and 3 both encode the same isoform (a).