

Product datasheet for **MC227438**

AA414768 (NM_001272033) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: AA414768 (NM_001272033) Mouse Untagged Clone
Tag: Tag Free
Symbol: AA414768
Synonyms: OTTMUSG00000017000
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227438 representing NM_001272033
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCGTCGGGGCTTAAGGCCGAGCTGGAGTTCCTAGCGTCCATCTTCGACCAGGACCACGAGCGGCTGC
GCATCGTAGCTGGCAGCTGGATGAGCTGCAGTGCAGTTCCTGGTGCCGCCGGCCCGCGGGCAGCCC
GCCCTCGCCGCGCTCACCTGCACTGCTCCATCACCGAGTCTACCCATCTTCGCCCGCATATGGTTC
GTGGACTCTGACGATCCAGATCTCACATCGATCCTGGAACGACTGGAAGATTGAAAGCCCAACAGTTCGC
TTCGTCAGCAGCTCAAGTGGCTGATCTGTGAACCTGCGAGCTTATAACCTGCCTGAGCACCTGGATAT
GGAGATGCTGGACCAGCCGGTCCCCATGGCTCAGAAGGCGACCATCGAAGAGGTGCCTTCAGAAGATGAA
GGTGAAGATGCAGCTGACGGGTGGAAGACTTGGACCAATTATGAGCTGAAAGAAGGAGGAGGAGGAGGAG
GAGGAGGAGGAGGAGAAGAAGAAGAAGAAGAGGCTAGTAGCGAGAAAAAGTCAGAGGATGAAGGGCT
TGAGAAAGAGCATTGGCAATATTAGAGAATATCAAGAAGAGCCAAAGGCAAGAGCATTGAAATGGTACC
GTGTCTGGCTCGGTGCATGCTTCAAACAGACTTATGAAAGAACTCAGAGATATATACAGATCGCAGAGCT
ACAAATCGGGGACTTATCCGTGGAACCTGATAAACGACAGCTTATGACTGGCAGCTTAAACTGAGGAA
GGTTGACCCCGATAGTCCCTTGTACGGTGACCTTCAGCTCTTGAATGAAAAGAAGGCATTGACTATATT
TTGCTTAACTTCTCTTCAAAGATAACTTTCCGTTTGATCCTCCATTTGTTGCGCTCGAGTTGCCTATTC
TCTCTGGAGGATGTTTTGGTGGAGGGCGTTGTGTATGAACTGCTCACAAACAGGGCTGGAGCAG
CGCCTACTCCATTGAATCAGTCATCATGCAGATCAATGCCACCTTAGTCAAAGGCAAGGCCAGGGTGC
TTTGGAGCAGATAAGAATCAATATAATCTAGAAACAGCCCAACAGTCTATGACTCCGTCGTCAGATGC
ATGAGAAAAAGGGCTGGTTCCTCCTCAAAGGAAGATGGG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001272033
Insert Size:	1164 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001272033.1</u> , <u>NP_001258962.1</u>
RefSeq Size:	1680 bp
RefSeq ORF:	1164 bp
Locus ID:	245350
Cytogenetics:	X A1.1