

Product datasheet for MC227502

Kcnj15 (NM_001271690) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnj15 (NM_001271690) Mouse Untagged Clone
Tag: Tag Free
Symbol: Kcnj15
Synonyms: 4930414N08Rik; AI182284; AI267127; IRKK; Kir4.2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227502 representing NM_001271690
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGTAGCCAGGTGGGAGAAGGGGAGCGAGGACGCGCCACTTACCCTGCAGAAGATTCTGACCTACAGA
 GTGGTCCAAAGATGGATGCCATTCACCTTGGCATGTCCAGTGCCCACTGGTGAAGCATACCAACGGGGT
 TGGACTCAAGGCCACAGACCCCGAGTCATGTCAAAGAGTGGGCACAGTAATGTGAGAATCGATAAGGTA
 GACGGAATCTATTTACTCTACCTCCAGGACTTGTGGACAACCGTCATCGACATGAAGTGGCGATACAAGC
 TCAACCCTATTTGCTGCCACCTTTGTGATGACCTGGTTTCTGTTTGGAGTGGTCTACTATGCCATAGCCTT
 TATTCATGGTGACTTACAACCTTGGGGAATCTAATTCCAACCACACACCCTGCATTATGAAAGTGGACTCT
 CTCACAGGAGCATTCTCTTTTCTTGGAACTCAGACAACCAATTGGCTACGGGGTCCGTTCCATCACAG
 AGGAGTGTCCCATGCTATCTTCTCTTAGTCGCCCAACTGGTCATCACCACATTGATTGAGATCTTCAT
 TACGGGGACCTTTCTGGCTAAAATTGCAAGACCCAAAAGCGAGCCGAGACCATTAAGTTCAGCCACTGT
 GCTGTATCAGCAAGCAGAATGGAAAGCTATGCCTGGTCATCCAGGTGGCCAACATGAGGAAGATCTCC
 TGATTCAGTGCCAGCTCTCTGAAAACCTCTGCAGACACAGTCAACAAAGAGGGAGAACGCATTCTCCT
 CAACCAGCCACTGTCAAATCCAGTGGACTCCTCTTCCGAGAGTCCCTTCTCATCTGCCCATGACC
 TTCTACCACGTGTTGGATGAGACAAGCCCCCTGCGGGACCTCACACCCAAAACCTAAAGGAGAAGGAGT
 TTGAGCTGGTGGTACTTCTCAACGCCACGGTGGAGTCTACCAGCGCCGTCTGCCAGAGCCGAACGTCTTA
 CATCCCGGAGGAGATCTACTGGGCTTTGAGTTTGTGCCTGTGGTTTCTCTCTCCAAAATGGAAAGTAT
 GTGGCTGATTTCAAGTCAATTTGAGCAGATCAGGAAGAGCCCGATTGTACCTTCTACTGTCCGATTCTG
 AGAAGCAGAAGCTTGAAGAACAGTACAGGCAAGAGGACCAGGGAGCGGGAGCTGAGGAGCCTCCTGCT
 ACAGCAGAGCAATGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001271690
Insert Size:	1209 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271690.1</u> , <u>NP_001258619.1</u>
RefSeq Size:	5138 bp
RefSeq ORF:	1209 bp
Locus ID:	16516
Cytogenetics:	16 55.86 cM
Gene Summary:	<p>Inward rectifier potassium channels are characterized by a greater tendency to allow potassium to flow into the cell rather than out of it. Their voltage dependence is regulated by the concentration of extracellular potassium; as external potassium is raised, the voltage range of the channel opening shifts to more positive voltages. The inward rectification is mainly due to the blockage of outward current by internal magnesium.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR compared to variant 1. Variants 1, 4, 6, 8, and 10 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>