

## Product datasheet for **MC227503**

### Kcnj15 (NM\_001271692) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kcnj15 (NM\_001271692) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Kcnj15  
**Synonyms:** 4930414N08Rik; AI182284; AI267127; IRKK; Kir4.2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227503 representing NM\_001271692  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGTAGCCAGGTGGGAGAAGGGGAGCGAGGACGCGCCACTTACCCTGCAGAAGATTCTGACCTACAGA  
 GTGGTCCAAAGATGGATGCCATTCACCTTGGCATGTCCAGTGCCCACTGGTGAAGCATACCAACGGGGT  
 TGGACTCAAGGCCACAGACCCCGAGTCATGTCAAAGAGTGGGCACAGTAATGTGAGAATCGATAAGGTA  
 GACGGAATCTATTTACTCTACCTCCAGGACTTGTGGACAACCGTCATCGACATGAAGTGGCGATACAAGC  
 TCAACCCTATTTGCTGCCACCTTTGTGATGACCTGGTTTCTGTTTGGAGTGGTCTACTATGCCATAGCCTT  
 TATTCATGGTGACTTACAACCTTGGGGAATCTAATTCCAACCACACACCCTGCATTATGAAAGTGGACTCT  
 CTCACAGGAGCATTCTCTTTTCTTGGAACTCAGACAACCAATTGGCTACGGGGTCCGTTCCATCACAG  
 AGGAGTGTCCCATGCTATCTTCTCTTAGTCGCCCAACTGGTCATCACCACATTGATTGAGATCTTCAT  
 TACGGGGACCTTTCTGGCTAAAATTGCAAGACCCAAAAGCGAGCCGAGACCATTAAGTTCAGCCACTGT  
 GCTGTATCAGCAAGCAGAATGGAAAGCTATGCCTGGTCATCCAGGTGGCCAACATGAGGAAGATCTCC  
 TGATTCAAGTCCAGCTCTCTGAAAACCTCTGCAGACACAGTCACCAAAGAGGGAGAACGCATTCTCCT  
 CAACCAGGCCACTGTCAAATCCAGTGGACTCCTCTTCCGAGAGTCCCTTCTCATCTGCCCATGACC  
 TTCTACCACGTGTTGGATGAGACAAGCCCTGCGGGACCTCACACCCAAAACCTAAAGGAGAAGGAGT  
 TTGAGCTGGTGGTACTTCTCAACGCCACGGTGGAGTCTACCAGCGCCGTCTGCCAGAGCCGAACGTCTTA  
 CATCCCGGAGGAGATCTACTGGGCTTTGAGTTTGTGCCTGTGGTTTCTCTCTCCAAAATGGAAAGTAT  
 GTGGCTGATTTCAAGTCAATTTGAGCAGATCAGGAAGAGCCGGATTGTACCTTCTACTGTCCGATTCTG  
 AGAAGCAGAAGCTTGAAGAACAGTACAGGCAAGAGGACCAGGGAGCGGGAGCTGAGGAGCCTCCTGCT  
 ACAGCAGAGCAATGCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001271692
<b>Insert Size:</b>	1209 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001271692.1</a></u> , <u><a href="#">NP_001258621.1</a></u>
<b>RefSeq Size:</b>	5106 bp
<b>RefSeq ORF:</b>	1209 bp
<b>Locus ID:</b>	16516
<b>Cytogenetics:</b>	16 55.86 cM
<b>Gene Summary:</b>	<p>Inward rectifier potassium channels are characterized by a greater tendency to allow potassium to flow into the cell rather than out of it. Their voltage dependence is regulated by the concentration of extracellular potassium; as external potassium is raised, the voltage range of the channel opening shifts to more positive voltages. The inward rectification is mainly due to the blockage of outward current by internal magnesium.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (8) differs in the 5' UTR compared to variant 1. Variants 1, 4, 6, 8, and 10 encode the same isoform (a). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>