

## Product datasheet for **MC227508**

### **AI987944 (NM\_001199330) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** AI987944 (NM\_001199330) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** AI987944  
**Synonyms:** MGC60751; MGC117553  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC227508 representing NM\_001199330  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGACCTATGATGATGTCCATGTAACCTCACTGAAGAAGAGTGGGCTTTGCTCAGTCCTTCCCAGAAGA  
 GGCTCTACAAGGATGTGATGCTGAAACCTACAAGAACCTCACTGCTATAGGTTACAATTGGGAAGACCA  
 TAATAGTGAACATTTTCAACATTCTAGAAGACATGGAAGGCATGAAAGAAGTCATACTGGAGAAAAACCT  
 TCTGTATATACTCAGTGTGGTAAAGTCTTTGCATATCGCAGTCATCCTCAGAATCATGATAGAATGCATA  
 CAGGAGAGAAAAGCCATGAAGGTATTCAACATGGTAAAGCCTTTGCACATCAAAGTCACCTTCAAAGGCA  
 TGAAAGAATTCATACCAGCAAGAAACCTTATGAATGTAATCAGTGTGGTAAAGCCTTTGCACATCGTAGT  
 TACCTTCGAAAGCATAAAAAGAAATTCATATTGAAGAGAAGCCCGATGAATGGAACCAATGTGGTCAAGGTA  
 ATCTCAAATACATAAAAGGACACATTCTGCACAGAAACCTATGAATGTAATCAGTGTGGAAAAGCCTT  
 CAGGCGGCACAATAGTCTCCAAGTACATAATAGAACACATACTGGAGAGAAAACCTACAAATGTGATCAA  
 TGTGATAAAGCCTTTACACAACAAGGTCATCTCGAATTCATAATAGAACACATACTGGAGAGAAAACCTT  
 ATGAATGCAATCAATGTGGTAAAGCCTTTGCAGTTCACAGTCTCCTCAAATACATAAAAGGACACATAC  
 TGCAGAGACACCCATGAATGTAATCAGTGTGGTAAAGCCTTTAGCGGCACAATACCCTCAAATACAT  
 AAAAGAACACACACTGGAGAAAAACCATACAAATGCGATCAATGTGATAAAGCCTTTTACAACAAGGTC  
 ACCTCAGAATACATAAAAGAACACATACTGGAGAGAAAACCTTATAAATGTAGTCAGTGTGGTAAAGCCTT  
 TGCACGCCATAGTCATCTTCTAATTCATGAAAGAATTCATACTGGAGAGAAAACCTTACAGATGTAATCAA  
 TGTAATAAAGCCTTTTACAATATGGTGGTCTTCAAATACATAAAAGCACACATACTGCTGAGAAGCCTT  
 ATGAATGTAATGAATGTGGTAAACTTTTTACAACACAGTTATCTCGAACACATAAAAGAACACATAC  
 TAGAGAGAAAACCTATAAA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001199330
<b>Insert Size:</b>	1212 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001199330.1</a></u> , <u><a href="#">NP_001186259.1</a></u>
<b>RefSeq Size:</b>	2948 bp
<b>RefSeq ORF:</b>	1212 bp
<b>Locus ID:</b>	233168
<b>Cytogenetics:</b>	7 B3