

Product datasheet for **MC227705**

Irx6 (NM_001301137) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Irx6 (NM_001301137) Mouse Untagged Clone
Tag: Tag Free
Symbol: Irx6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC227705 representing NM_001301137
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTTTTCACCTTTGGACACCCGTATGGCAGCACATCCCAGTTTCTGGTGTCTGCAAGTTCAGTG
 CCACTTGTGCGAAACGGCCCTCGGCCAGTGTGAGATGTGGCCTCAGCCTCCACCTCTGCCTCTACTCT
 TTGCTGTACACCCATGACAGTCGGCTGCTGGCAGTGTGCGCCAGAGCTGGGAGCTGCCTTGGGCATC
 TATGGAGCACCCATGACAGTCGCCAGAGTTACCCTGGGTACCTGACCTATGGCCAGAGCCACCCACAC
 TGTGTGGCGCACTGAATCCTCAGTATGAGTTCAAGGATGCTGCAGGAAGCTTTGCCCCAGCCTGACCCA
 GCCAGGGGCTACTACCCCTATGAAACAACCTCTGGGGCAGTATCAATACGACAGGTATGGCGGAGTGAA
 TTGAGCAGTGTGGCCGAGGAAGAATGCCACAAGGGAGAGCACCAGCGCACTGAAGGCCTGGCTCCACG
 AGCACCGCAAGAACCCGTACCCACCAAGGGCGAGAAGATCATGCTGGCCATCATCACCAAGATGACCC
 CACCCAGGTGTCCACCTGGTTCGCCAACGCGCCGCGCCCTCAAGAAGGAGAACAAGATGACATGGGCG
 CCCAAGAACAAGGGAGGGGAGGAGAGGAAAGCAGACAGTGGAGAAGACTCTCTGGGATGCCTGAATGGT
 ACACCAAAGATGCTACTGCCAGCCAGGAAGCCCGGGGCTGAGGCTGAGTGACCTGGAAGACCTGGAGGA
 AGAGGAGGAAGAGGAAGAAGCAGAGGAGGAGGCAGCGGTCTCAGCAGCACGAGGCTGGCGGATTTTCAG
 AAGAGCACGCAACCCCTGCCTGCTCCCTGCGCCGCGCCAGGAGAGATGCTTGGAGAGCAGAGAGTGCG
 GCCTGGGGCTACCCCGTTTCTCTTTTACTGAGGCCCAACAATCAGGGGAAGCTGACTTATCACAGCAGA
 GCCAGGGCGCCCCACGATGATCTTACACTACCCAAGTGGCCACAAACCCCGCATTGGTCCCTTAGCTCAC
 ACTGCGGCAGCCAGCGTGTGAAAAGTGTCCCTCAACTCCGCCAGGGCACAAGTCCAGAGTGCCACA
 TGATTCAGACAGCCAGTCTATCAGGCGACTCCTGGTACCCAGAGACTCCGAAGGCGAAGAGGATTC
 TCCTGCAGCCAAAGCTTTTGGAACTCCACGTTACCCCTGCAGGGGCTGCCACTGAATGTGCACCATAC
 CCGAGGCGGAGAGAGCCTGAAGTGGGTTCCAGTACCCATCTGGAGCAGAAGG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001301137
Insert Size:	1317 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001301137.1</u> , <u>NP_001288066.1</u>
RefSeq Size:	3366 bp
RefSeq ORF:	1317 bp
Locus ID:	64379
UniProt ID:	<u>Q9ER75</u>
Cytogenetics:	8 44.99 cM