

## Product datasheet for **MC227830**

### 2410141K09Rik (NM\_001290196) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	2410141K09Rik (NM_001290196) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	2410141K09Rik
Synonyms:	4930482J15Rik; Ayu21-35; Gt(pU21)35Imeg
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**Fully Sequenced ORF:** >MC227830 representing NM\_001290196  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGATGCAATGACTTTTGACGATGTGCATGTGAACTTCACTAAAGAAGAGTGGAAATTTGCTGGAACCTT  
CCCAGATGAATCTCTACAAGACGTGATGTTGGAGACCTACTGGAACCTTACTTCTATAGGCTACAATG  
GGAAGACCATCATATTGAAGAACCATGTCAAAGTTCTCGAAGACATACAAGGCATGTCAGGAGTCATAGT  
TGGAAGAAACTCTATGAATGTAATCAATGTGGTAAAGCCTTTACCTCTCACAGTAGTCTCAGATATCATA  
AAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCATTGTGGTAAAGCCTTTGCAAGACCTGCTCA  
TCTCCAAAGACATAAAAGAACCATACTGGAGAGAAACCTTATGAGTGAATCAATGTAGTAAAGCCTTT  
TCCTGTACAGTGGTCTCAGATATCATAAAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCAAT  
GTGGGAAAGCCTTTGCAAGACCCGCTGATCTCCAATGTCATAAAAGGACCATACTGGAGAGAAACCTTA  
TGATGTAATCAATGTGGAAAGCCTTTGCAAGACCCGCTAATCTCCAATGTCATAAAAGAACCATACT  
GGAGAGAAACCTTATGAATGTAATCAATGTGGTAAAGCCTTTCTTGTACAGGTTGGTCTCAGATATCATA  
AAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCAATGTGGGAAAGCCTTTGCAAGACCCAGTCA  
TCTCCAATGTCATAAAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCAATGTGGTAAAGCCTTT  
TCCTCTCACATTAGTCTCAGATATCATAAAAGAACCATACTGGAGAGAAAGCCTTATGAATGTAATCAAT  
GTGGTAAAGCCTTTCCGGTCACAGTGGTCTCAGATATCATAAAAGAACCATACTGGAGAGAAACCTTA  
TGAATGTAATCAGTGTGGTAAAGCCTTTGCAACACCCAGTCATCTCCAATGTCATAAAAGAACCATACT  
GGAGAGAAACCATATGAATGTAATCAATGTAGTAAAGCCTTTGCAAGACCCGCTCATCTCCAACATCATA  
AAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCAATGTGGGAAAGCCTTTGCAAGGCTCAGTCA  
TCTCCAATGTCATAAAAGAACCACACCCGGAGAGAAACCTTATGAATGTAATCAATGTGGGAAAGCCTTT  
TCCTGTACAATAGTCTCAGATATCATAAAAGAACCATACTGGAGAGAAACCTTATGAATGTAATCAGT  
GTGGTAAAGCCTTGCTACATCCACTCATCTCCAATGTCATAAAAGAACCATACTGGAGAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM\_001290196

**Insert Size:** 1395 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001290196.1](#), [NP\\_001277125.1](#)

**RefSeq Size:** 1491 bp

**RefSeq ORF:** 1395 bp

**Locus ID:** 76803

**Cytogenetics:** 13 B3