

## Product datasheet for **MC228019**

### Zfp708 (NM\_001292066) Mouse Untagged Clone

#### Product data:

|                    |  |
|--------------------|--|
| Product Type:      | Expression Plasmids                        |
| Product Name:      | Zfp708 (NM_001292066) Mouse Untagged Clone |
| Tag:               | Tag Free                                   |
| Symbol:            | Zfp708                                     |
| Synonyms:          | Rslcan11                                   |
| Vector:            | pCMV6-Entry (PS100001)                     |
| E. coli Selection: | Kanamycin (25 ug/mL)                       |
| Cell Selection:    | Neomycin                                   |



[View online »](#)

**Fully Sequenced ORF:** >MC228019 representing NM\_001292066  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTTGGAGAATTACAGCAACCTTGTGTTCTGGTCTTCTCTAAGCCATACCTGGTCACATTTCT  
TGGAGCAGAGTTCAGATCTGTGGAGTGTGAAGAGGAAGGCATCAGTGGCCATCCATCTAGGATCAAAGC  
CCATAAATGTAAGAATGTGGGAATCCTTTGATAGGAAGTCACTCCTTATTCAACATCAAAGAATTCAT  
ACTGGAGAGAGGCCCTACAAATGTGAAGAATGTGGCATGGCCTTCACCGGTCGCTGCATTCTTTCTAAGC  
ACCAGCGAATCCATACTGGAGAAAAACCTACAAATGTAAGAATGTGGCAAAGCTTTCAGCTGTTCTTC  
TAACCTTCACCAACACCAACAAATTCATAGGGAAGAGAAACCTATAAATGCGATGACTGTGGGAAGGCC  
TTTAGCTGTTCTCATATCTGTATAAGCATCAGCGAATCCATACTGGCATGAAACCCTACAAATGCAAGG  
AATGTGGCAAGGCCTTTAACTGTTCTGCAACCTTATTTACCACCAGAGAATTCATACTGGAGAGAAACC  
ATATAAGTGTGGTGAATGTGGGAAGGCCTTACCATTTATTCAACATTTATGAATCACCAGCGAACACAT  
AGTGGAGAAAAACCGTATGAATGCAAGACTGTGGGAAAGCCTTACTAATTGTTATAATCTAATTC AAC  
ACCAAAGAATTCATACTGGAGAGAAACCTATGCATGTGAAGAATGTGGCAAATCTTTAACTATTCTTC  
CCATCTTAAACAACAAAAAATCCATACTGGAGAGAAACCATATAAATGTGAAGTGTGTGGCAAAGCT  
TTAACTGTTCCGATTATCTTGTAAAGCATCAAAGAATTCATACTGGAGAGAAACCTACAAATGTGAAG  
TATGTGGCAAAGCTTTTAGCTTTTCTACATATCTGCATAAACATCAAAGAATCCACACTGGAGAGAAACG  
CTACAGATGTGAAGAATGTGTAAGCTTTTACTAATTATTACAGGCTTATTGTTACCGAAGAGTACAT  
ACTGGAGAGAAACCTACAAATGTGAAGAATGTGGCAAGGCCTTAGTGTTCACACATCACTCTCTAAAC  
ATCAGAGAATTCATACGAGAGAGCAACCTTACAAATGTGATGAATGTGGAAAGACCTTTAGTGTTCAC  
AACATTTCTAAACATCAGATAATGCATACTAGAGAGAAACCTACAAATGCAAAACACTGCGACAAGACC  
TTCAAATATCCTTACCCTTGTCAACATGAAAGAATTCACACTGGAGAGAGACCCTACAAATGTGAAG  
TTTGTGGCAAGGCCTTTACTTCTTCTTCAAATCTTAAATATCATTGGAGACTTCATACTGGCGAGAAGCC  
CTATAAATGTGAACAATGTGGTAAAGCCTTTAAAAATAGTTCAAATGTCCTCGACACTACAAAATTCAT  
ACTAGAGAAAATCCAGAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001292066

**Insert Size:** 1491 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001292066.1](#), [NP\\_001278995.1](#)

**RefSeq Size:** 2436 bp

**RefSeq ORF:** 1491 bp

**Locus ID:** 432769

**Cytogenetics:** 13 B3