

Product datasheet for **MC228040**

Btbd35f13 (NM_001282034) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Btbd35f13 (NM_001282034) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Btbd35f13
Synonyms:	100040722; Gm2927
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC228040 representing NM_001282034
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGGCTTTTAGTCAGCAGGTTCTTGAGATGCAGGGATTCCAGTCTGTTAGAGCCACAGCCAGAAGCCA
 TAGCCGGAGCCAGCTACATTCTGGCAGTCGCAAGCGAAAAAGAAACAGTTTGGAGGAGTTGGCAACAAG
 TTCTAATGTTTCATGGACCTCAAAACCAGGGAATGTATCCACATCAAGTTCTCAACTACATCTACTGGAAA
 AGGGTTAAGATCTCATCTAATGATGCTTATCAAAACTATTTTTGGACGGGCATGATAGCGACATTA
 TCCGTGCTCTGGAAAAACATGGTGTACACAAAGTATTTTTATGTCAGTCAGGCTACTTTGCTAACAT
 TCTCAAAGGTAAGTGGAGAGAATCACACCATGGTGTATAAATCTGATCATTAAAGATGAGGATATTGAT
 ACCCGATCTGCAATTTGTGTTGGTGTGTTGTACACAGATGCGGATTTGTCAATAACACCTCTGGAAG
 TTCTCAAGTTTGGCAGCAGCATGCCTGCTTCGGGTGGATCGAGTAATTCAGCAGTGTGAAGGAATCAT
 GAAAGAACTATCAACAGGAACACTGTGTCTCTATTATTTGGCAGCAGAAACCTATAGATTA
 AAGCTGTAAGAGCAGATGCTTTGAATGGCTTCTTTGCAATTTGATGGTACATCCAAGTGTGGC
 ACTTTACAAGG AAGTAGACTTGAAGTTGATGTATGTTCTAGCACTGTCTTCTGACTTACTAGTCATGCA
 AAAAGGAGATTGATGTATATAACCACACTAAAAATATGGATGTTCTTTATCTTAATCCATGCTGGA
 AACGGAACCATGAAACAGCTTTTACAACACGCAAACTGGCTTTCCACCCACATGGC
 ATATGTTGATAACATCAGTTTTCTTGAAA GTGAAGAAGGACTAATATTTCAACCAAGTGT
 TAAAAAGCTGAGATTTTCAGCACATTATCTGTGACTTGAC TTCCACAACTATCTTGAACAAGAT
 CGACTAATACCTATGGCATGGTTGTCACCCATTTACAACAACAG TGGTTGACTTTGCTGCGAACACA
 AAGAATATGGGGTAATGGACCACAAGTATCAATGAACAAGAAGTTG AAGAATGCACCATGAGGTGT
 GGTACAATGATCCCAAGGATGGAAGATATACTTGAAGTGGTCAAGTTGG ACGACTTGGCTTTCC
 TTTACGTGTGACCTTTACCAGGCAGTGTGTAATTTAAGGCAACGGTGTCAAGGG TGTGATG
 GTTCTGCTTGCACAACCATATCCGAAATGTCATATTCAGAATAACTTTGGTGTGTTTTGATT
 CCAACAAAAGAGTAACTTTCAGAAAGACAACAGGTTATAAAATCCTCACCTTTGAATATAAGGAGG
 AGCA AATTGTAATGAAATTGGATAGTGATGTTCTAACCTTCCCTATGTGTATATTCTGCAATTTCT
 TTTTGTA AACCTAGGAAATGCAGAAAAACAAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_001282034
- Insert Size:** 1497 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001282034.1](#), [NP_001268963.1](#)

RefSeq Size: 1943 bp

RefSeq ORF: 1497 bp

Locus ID: 100040722

Cytogenetics: X