

Product datasheet for **MC228163**

Cyp4f14 (NM_001204333) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyp4f14 (NM_001204333) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyp4f14
Synonyms:	1300014O15Rik; cyp4f-14
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC228163 representing NM_001204333
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCACAGCTGAGCCTGTCCTGGTTGGGACTGGGGCCCGAGGTGGCCTCCCATGGAAGACCCTGCTAC
 TGCTTGGGGCCTCCTGGATCCTGGCCCGATTCTGATCCAGATCTATGCTGCCTATAGGAACTATCGTCA
 CCTCCATGGCTTCCCTCAGCCCCCAAACGGAAGTGGCTTATGGGTCACGTGGGCATGGTAACTCCTACA
 GAGCAGGGCTTGAAGGAATTGACTCGTTTTGGTGGGTACCTACCCCAAGGTTTTTTGATGTGGATTGGAC
 CAATGGTCCCTGTTATTACTCTATGCCATTCTGACATCGTCCGATCTATCCTCAATGCCTCAGCCGCTGT
 TGCCTCAAGGATGTGATATTCTACAGCATCTAAAGCCTTGGCTGGGGGATGGGCTCTTGGTGAAGTGT
 GGTGACAAGTGGAGCCGCATCGTCGCATGCTGACACCCGCTTCCATTTCAACATCCTGAAGCCCTATG
 TGAAGATTTTCAACGATAGCACCAACATCATGCACGCCAAGTGGCAACGGTTGATTTAGATGGCAGTGC
 CCGTCTGGACATGTTTGGACAGCTCAGCCTGATGACTCTGGATAGTCTGCAGAAATGTGTCTTCAGCTTT
 GACAGCAACTGTCAGGAGAAGTCCAGTGAATATATTGCTGCCATCCTGGAGCTCAGTGCCTAGTGGCCA
 AAAGGCACCAGCAGCCTTTGATGTTTCATGGACCTGCTGTATAATCTCACCCCTGATGGTATGCGCTTCCG
 TAAGGCCTGCAACGTGGTGCATGAATCACTGATGCTGTATCCGGGAGCGGCACCCGACCCCTCCAGAT
 CAGGGTTTGGATGACTTCCCTAAGTCTAAGGCCAAGTCCAAGACTTTGGATTTTCATTGATGTGCTGCTGC
 TGAGCAAGGATGAAGATGGAAAGGAGCTGTGAGATGAGGACATCCGAGCAGAGGCTGACACCTTCATGTT
 TGAGGGCCATGACACGACAGCCAGTGGACTCTCCTGGATCCTCTACAACCTGGCAAGGCACCCAGAATAC
 CAGGAGCGCTGCCGCGAGGAAGTTCAGGAGCTGCTGAGGGCCGAGAGCCTGAGGAGATTGAATGGGACG
 ACCTGGCCAGCTACCTTCTGACCATGTGCATCAAGGAGAGTCTGCGACTGCATCCCCAGTACCGGT
 CATCTCCCGATGCTGCACCCAGGACATTCTGCTCCAGATGGCCGACCATTCCCAAAGGTATCATCTGT
 CTCATCAGCATCTTTGGAATTCATCACAACCCATCCGTGTGGCCAGACCCCTGAGGTCTATGACCCCTTCC
 GCTTTGACCCTGAAAACATCAAAGACAGCTCACCTCTGGCATTATTCCTTCTCGGCGGGACCCAGGAA
 CTGCATAGGACAGACTTTCGCCATGAGTGAAGTGGCACTGGCGCTGACGCTGCTGCGCTTCCGA
 CTCCTACCGGATGACAAGGAGCCGCGCAGGCAGCCAGAGCTGATCCTGCGCGCAGAGGGCGGGCTGTGGC
 TGGGGTGGAGCCGCTGAGCGCGGGCGCTCAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul

ACCN: NM_001204333

Insert Size: 1575 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204333.1](#), [NP_001191262.1](#)

RefSeq Size: 2054 bp

RefSeq ORF: 1575 bp

Locus ID: 64385

UniProt ID: [Q9EP75](#)

Cytogenetics: 17 B1

Gene Summary: Cytochromes P450 are a group of heme-thiolate monooxygenases. Catalyzes the omega-hydroxylation of LTB4 (By similarity).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, 4 and 5 encode the same protein.