

## Product datasheet for **MC228334**

### Arntl2 (NM\_001289679) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Arntl2 (NM_001289679) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Arntl2
Synonyms:	4632430A05Rik; bHLHe6; BMAL2; CLIF; MOP9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >MC228334 representing NM\_001289679  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGAGTTTCCAAGGAAACGCAGAGGCAGAGATCCCAGCCACTCCAGTCAGAATTCATGACAGACACAA  
 CAGTGGAAAGTCTTCCCCAGAATCCCTTTGCCTCTCTTTCAACAAGAACAGGAGTATCAGCGCCAG  
 TGGCATCAGGGAAGCTCACAGCCAGATGGAAAAGCGTCGGAGAGACAAGATGAACCATCTGATTCAGAAA  
 CTGTCTATGATCCCTCCACACATCCCCACGGCCACAAAAGTGGACAAGCTCAGCGTCTTGAGGAGGG  
 CGGTGCAGTACTTGAGGTCTCTGAGAGGCATGACAGAGCTTTACTTAGGAGAAAAGTCTAAACCTTCATT  
 TATTCAGGATAAGGAACTCAGTCACTTAATCCTCAAGGCAGCAGAAGGCTTCTGTTTGTGGTTGGATGC  
 GAAAGAGGGAGAATTTTTACGTTTCTAAGTCTGTCTCCAAAACACTGCGTTATGATCAGGCTAGCTTGA  
 TAGGACAGAATTTGTTTACTTCTACACCCAAAAGACGTCGCCAAAGTAAAGGAACAATTTCTTGTA  
 TGGTTCACCAAGAGAGAAACCTATAGACACAAAACCTCTCAGGTTTACAGTCACCCCTACACTGGGCGA  
 CCACGCATGCATTCTGGCTCCAGACGATCTTTCTTTAGAATGAAGAGCTGTACCGTCCCTGTCAAAG  
 AAGAGCAGCCATGCTCGTCTGCTCAAAGAAGAAAGACCATAGAAAATCCACACCGTCCATTGCACTGG  
 ATACTTGAGAAGCTGGCCTCTGAATGTTGTTGGCATGGAGAAAGAGTCCGGTGGTGGGAAGGACAGCGGT  
 CCTTTACCTGCCTGTGGCTATGGGACGGTGCATCCATACATTGTCCCTCAAAGAGTGGCAAGATCA  
 ACGTGAGACCGGTGAGTTCATAACTCGCTTCGCAATGAACGGGAAATTCGTCTATGTTGACCAAAGGGC  
 AACGGCAATTTTAGGATACCTGCCTCAGGAACTTTGGGAACTTCATGTTATGAATTTTCATCAGGAT  
 GACCACAGTAGTTGACTGACAAGCACAAGCAGACTCATACAAATTCAGAGTGAAGGATGGTGCCTTCG  
 TGACTCTGAAGAGTGAGTGGTTCAGCTTCACAAACCCTTGACCAAAGAGCTGGAGTACATTGTGTCTGT  
 CAACACATTGGTTTTGGGGCGCAGTGAGACCAGGCTGTCTTTGCTTCATTGCGGGCGCAGCAGCCAGTCC  
 TCCGAAGACTCATTTAGACAATCCTGCATCAATGTGCCCGGTGTATCCACGGGGACCGTCTTGGTGCTG  
 GGAGTATTGGAACAGATATTGCAAATGAGGTTCTGAGTTTACAGAGATTACACTTTCATCCCCAGAAGA  
 TGCAAGCCCTTCAGAAGAAGTGAGAGATGACTGCAGTGTAAATGGTGGGAATGCCTATGGCCTGCATCC  
 ACTAGGGAGCCTTTTGCAGTGAGCCCTTCTGAAACAGAGGTCTGGAGGCTGCCAGGCAACACCAGAGCA  
 CTGAACCCGCCACCCTCACGGACCACTCCCGGTGACAGTGCCAGCTGGGTTTTGATGTCCTGTGTGA  
 CAGTGACAGCATAGACATGGCTGCATTCATGAATTACCTCGAAGCAGAGGGGGCCTGGGTGACCCTGGG  
 GACTTCAGTGACATCCAGTGGGCACT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** SgfI-MluI
- ACCN:** NM\_001289679
- Insert Size:** 1710 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001289679.1, NP_001276608.1</u>
<b>RefSeq Size:</b>	2593 bp
<b>RefSeq ORF:</b>	1710 bp
<b>Locus ID:</b>	272322
<b>UniProt ID:</b>	<u>Q2VPD4</u>
<b>Cytogenetics:</b>	6 G3
<b>Gene Summary:</b>	<p>This gene encodes a basic helix-loop-helix transcription factor belonging to the PAS (Per, Arnt, Sim) superfamily. The PAS proteins play important roles in adaptation to low atmospheric and cellular oxygen levels, exposure to certain environmental pollutants, and diurnal oscillations in light and temperature. This protein forms a transcriptionally active heterodimer with the circadian Clock protein, the structurally related Mop4, and hypoxia-inducible factors, such as Hif1alpha. Consistent with its role as a biologically relevant partner of circadian and hypoxia factors, this protein is coexpressed in regions of the brain such as the thalamus, hypothalamus, and amygdala. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Feb 2014]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site, compared to variant 1. The encoded isoform (2) is shorter than isoform 1.</p>