

## Product datasheet for **MC228421**

### **Rtl5 (NM\_001278534) Mouse Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Rtl5 (NM_001278534) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rtl5
Synonyms:	6430402L03Rik; Mar5; Mart5; mKIAA2001
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001278534.1](#), [NP\\_001265463.1](#)

**RefSeq Size:** 4330 bp

**RefSeq ORF:** 1800 bp

**Locus ID:** 331474

**UniProt ID:** [Q5DTT4](#)

**Cytogenetics:** X D