

Product datasheet for **MC228439**

Papss1 (NM_001289478) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Papss1 (NM_001289478) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Papss1
Synonyms:	AI325286; Asapk; SK1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC228439 representing NM_001289478
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGCAGAGAGCAACCAACGTACCTATCAAGCCCACCATGTCAGCAGGAATAAGAGAGGGCAGGTGGTGG
 GGACCAGAGGTGGCTTCCGTGGTTGCACAGTTTGGCTGACAGGATTGTCCGGAGCAGGGAAGACGACCGT
 GAGCATGGCTCTGGAGGAGTACCTGGTGTGCCACGGCATTCCGTGCTACACTTTGGATGGTGACAACATC
 CGCCAAGGACTCAATAAGAACCTCGGCTTCAGTCTGAGGACAGAGAAGAGAACGTTCCGCCGATAGCTG
 AGGTGGCGAAGCTGTTTGCAGATGCTGGCTTAGTGTATCACCAGCTTTATATCGCCTTACACACAGGA
 TCGCAACAACGCAAGGCAGATTCATGAGGGTGAAGCTTGCCCTCTTTGAAGTTTTTGTGATGCTCCT
 CTGCATGTCTGTGAGCAGAGGGATGTCAAAGCCCTCTACAAGAAGGCGCGGCAGGGGAGATAAAAGGCT
 TCACTGGCATCGATTCTGAGTATGAGAAACCCGAGGCCCGGAGCTGGTGTGAAAACGGATTCTGTGA
 CGTCAACGACTGGTCCAGCAGGTTGTGGAGCTTCTTCAGGAACGGGACATCGTCCCTGTGGATGCTTCC
 TATGAAGTGAAAGAGCTATATGTGCCAGAGAATAAACTTCACTTGCCAAAACCTGATGCAGAAGCCTTAC
 CAGCCCTGAAAATCAATAAAGTGGATATGCAGTGGGTGCAGGTTTTGGCAGAAGTTGGGCGACTCCTCT
 GAACGGCTTATGAGAGAGAGGGAGTACTTGCAGTGCCTCCATTTGATTGTCTTCTGGATGGAGGCGTC
 ATCAACTTATCGGTGCCTATAGTTCTGACAGCTACGCACGAGGATAAGGAGAGGCTGGACGGCTGCACCG
 CGTTTCGCTCTGGTGTATGAGGGCCCGCGTGGCCATCCTTCGGAATCCTGAATTTTTGAGCACCAGGAA
 AGAGGAGCGGTGTGCCAGACAGTGGGGAACAACATGCAAGAACCACCCCTACATCAAGATGGTTCTGGAA
 CAAGGGGATTGGCTGATTGGAGGAGATCTCAAGTCTGGACCGGATTTACTGGAATGATGGTCTTGATC
 AGTACCGCCTTACCCCGACGGAGCTCAAGCAGAAGTTAAAGATATGAACGCTGATGCTGTCTTTGCATT
 TCAGTTGCGCAACCCAGTGCACAACGGCACGCTCTGTTAATGCAGGATACCCACAAGCAGCTTCTGGAG
 AGGGGCTACCGGCGCCCTGTCCTGCTCCTCATCTCTTGGTGGCTGGACGAAGGATGACGATGTCCCTC
 TGATGTGGCGTATGAAGCAGCACGCTGCAGTGTGGAGGAGGCATCCTGGATCCTGAAACGACAGTGGT
 GGCCATTTTCCCGTCTCCTATGATGTATGCTGGGCCAACCGAGGTCAGTGGCACTGCAGGGCGCGGATG
 GTGGCCGAGCCAATTTTATATTGTTGGACGAGACCCTGCTGGCATGCCTATCCGGAGACAGGGAAGG
 ACCTCTATGAGCCAACACATGGTGCCAAAGTGTGACGATGGCCCCAGGCCTGATTACCTTGAAATCGT
 TCCCTCCGAGTTGCAGCTTACAACAAGAAAAAGAAGCGGATGGACTACTATGACTCTGAGCACCAGGAA
 GACTTCGAGTTTATTTTCAGGAACGAGGATGCGCAAGCTGGCACGGAAGGCCAGAAACCTCCAGAGGGCT
 TCATGGCCCCCAAGGCCTGGACTGTGCTGGTAGAGTACTACAAGTCTTAGAGAAAAGCCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** SgfI-MluI
- ACCN:** NM_001289478
- Insert Size:** 1812 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001289478.1](#), [NP_001276407.1](#)

RefSeq Size: 2794 bp

RefSeq ORF: 1812 bp

Locus ID: 23971

Cytogenetics: 3 61.05 cM

Gene Summary: Bifunctional enzyme with both ATP sulfurylase and APS kinase activity, which mediates two steps in the sulfate activation pathway. The first step is the transfer of a sulfate group to ATP to yield adenosine 5'-phosphosulfate (APS), and the second step is the transfer of a phosphate group from ATP to APS yielding 3'-phosphoadenylylsulfate (PAPS: activated sulfate donor used by sulfotransferase). In mammals, PAPS is the sole source of sulfate; APS appears to be only an intermediate in the sulfate-activation pathway (PubMed:7493984). Required for normal biosynthesis of sulfated L-selectin ligands in endothelial cells (By similarity).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (3) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (b) is shorter at the N-terminus compared to isoform a. Variants 2 and 3 both encode the same isoform (b).