

Product datasheet for **MC228710**

Adgrf4 (NM_001289501) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adgrf4 (NM_001289501) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Adgrf4
Synonyms:	4632435A09Rik; AV239010; Gpr115
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC228710 representing NM_001289501
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAGCCCTGGATTGCCATGGTTTGTGTTTTAGTATTTTTCTGACAACAGAATGTTCTCATTCCAAC
 CCAAGACTCACAGAAAAGATGAAGATAAATTTCAAATATCTCTACAGAAGCATGAGTTTAGACCAAGACA
 GGGTAAATGTGATGGACTGTGCAGCTCTTCTTCCAGCTGTAACCCAGAGCTGTCCTTGGAACCTCCCGTGA
 GAAATCGTATTCACATGTAATCAAAACAATGGCAAAAAACAATCGAAACATGCACAAGTCTCTGTGG
 ACACCTTGTTCAGCGCATACATCCTGCAGCTAGCCTATCCTTGGCGTCTTCTTGTTTTTCCCATGAG
 CTTAATTGGCAATGCAGCACCAGTCCATATTGAAATGTATTTCAAGGAATCCAAAAGTACTGCCCGAG
 GACTACGTTTGCATTGTTGATGCTGTGAAATCCTCTGCAGTTACATCTGGAAATATTGCATTATAGTAG
 AGTTATTAATAAATATCTCTTCAAACCTGCAAACCTTCTGGAATACATGATAATGTGAATTGGAAGAAAAT
 GAAGAATATGGCAAAGTGGCAACCATATCCTTGGCCCCACTGCCATTTCAAATGGGCTTTCATTGCC
 AACAAAAATGCCAGCTCAGATCTGTTGGAGTCGGTGAACCTTTTTGCCAAGAAATCCAAATCCAAGGTA
 AATCTGAGAGCATCGTAGATGAACCTTTTATTCAGACGAAAGGCTCTCGTATCTCTCACAGTTCCCTCAGA
 GCACAGCCTCAGCCTCTCTGTACCCAGGTACAATGCTACAGAAGATGCTTAGTGGTGATCGAAATTCCT
 AGGCAAGCATTACAGGAGCTGTCGTTCAATGCCTCCCAAGCCATTGTTGTCGCTTTTCCAACCTTGGGGG
 CCATTCTAAAAGAAGTCCATCGGCCAAATACAAATCTCCAAAAGCCTATAGATGACCTCATCCTATCACT
 GGTTTTGCCAGAAGGTCTGAACGAAATCATACTCACCTTTGATAAAATCAATAAATCCCAGAGCACCAGC
 AGCCAGTGTGTCAGCTGGGATCCTGCAACAGGGCAGTGGGATGAGAGCCCATGTACAGTGATGTCAGACA
 TCAATAGCACAGTGAAATGCCGGTGCAGGCACACCAAAGCTGTGACGTCCTTCTCCATCCTCATGTCTTC
 TAAGCCTGTAAAGAATACGATTCTGAACCATATCACCTTCATTGGACTCAGCATCTCCATCTTTAGCTTG
 GTGCTCTGCCTGGTCATTGAGGCCATCGTCTGGTCCCAGGTTGGTTTACAGAGATATCCTACATGCGCC
 ACGTGTGCATCGTGAACATTGCAGTGTCCCTCCTGACTGCCAATGTGTGGTTCATCATTGGCTCCAACCT
 CAGTGCCAACGTTCCAGGAGGATCACAAGTGGTGTGGCTGTGACATTTCTGCCACTTCTTCTCCTC
 TCCCTGTTCTTCTGGATGCTCTTCAAAGCGTGTCTCATAGTCTATGGAATACTGGTTGTCTTCCGAGGA
 TGATGAAGTCTCGGATGATGGCCATTGGCTTTGCCATTGGCTATGGGTGTCCCCTGTGATCGCGGTCAT
 CACTGTTACTGTCACAGAGCCTGGGAAGGCTACACCAGAAAAGATGCCTGTTGGCTTAATTGGAACCAA
 ACCAAAGCCCTTTTTGCTTTTGCTATCCCAGCCTTGGCCATTGTGGCTGTGAATCTCCTTGTGGTCTGG
 CAGTTGCTATCAACACACAGAGGCCCTTGATTGGCAGTTCAAAATCTCAGGATATGGCCATTGTTTTAG
 GATCAGCAAAAATGTTGCCATCCTACCCCACTTCTGGGACTGACCTGGGGTTTTGGTCTAAACCCCTG
 CTAGAAGGCGTACACCTGGTATTCATATAATCTTCGCTTGTCAATGCTTCCAGGGCTTTTTCATCT
 TGCTGTTTGGTACCATTATGGATCACAAGATAAGAGATGCTCTGAGGATGAGGGTATCTTCATTGAAGGG
 GAAATCAAGGGCAGCAGAGAAGGTGTCATTAAGTCCAGCCAATGGATCCAGAATACTGAATCGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul

ACCN: NM_001289501

Insert Size: 2097 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001289501.1</u> , <u>NP_001276430.1</u>
RefSeq Size:	3021 bp
RefSeq ORF:	2097 bp
Locus ID:	78249
UniProt ID:	<u>Q9D2L6</u>
Cytogenetics:	17 B3
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, and 4 encode the same protein.