

Product datasheet for **MG200028**

Polr2k (BC028543) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Polr2k (BC028543) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Polr2k
Synonyms: MafY, RPB12, RPABC4, RPB7.0, RPB10alpha
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200028 representing BC028543
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACGCTCAGAAAGATGTTCAACCACCAAAGCAGCAGCCAATGATATATATTTGTGGAGAGTGTCA
 CCGAAAATGAGATAAAGTCCAGGGATCCCATCAGATGCAGAGAATGTGGATACAGAATAATGTACAAGAA
 AAGGACTAAAAGATTGGTGGTTTTGATGCTCGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200028 representing BC028543
 Red=Cloning site Green=Tags(s)

MDAQKDVQPPKQPMIYICGECHEIKSRDPIRCRECGYRIMYKKRTRKRLVVFAR

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC028543

ORF Size: 176 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC028543 , AAH28543
RefSeq Size:	480 bp
RefSeq ORF:	176 bp
Locus ID:	17749
Cytogenetics:	15 B3.1
Gene Summary:	DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Common component of RNA polymerases I, II and III which synthesize ribosomal RNA precursors, mRNA precursors and many functional non-coding RNAs, and a small RNAs, such as 5S rRNA and tRNAs, respectively (By similarity). [UniProtKB/Swiss-Prot Function]