

Product datasheet for **MG200066**

Defb1 (NM_007843) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Defb1 (NM_007843) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Defb1
Synonyms: AW260221; BD-1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200066 representing NM_007843
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAACTCATTACTTTCTCCTGGTGATGATATGTTTTCTTTCTCCAGATGGAGCCAGGTGTTGGCA
TTCTCACAAGTCTTGACGAAGAACAGATCAATACAAATGCCTTCAACATGGAGGATTCTGTCTCCGCTC
CAGCTGCCATCTAATACCAAACACAGGGAACCTGTAAACCAGATAAGCCCAACTGTTGTAAGAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200066 representing NM_007843
Red=Cloning site Green=Tags(s)

MKTHYFLLVMICFLFSQMEPGVILTSLGRRTDQYKCLQHGGFCLRSSCPSTNKLQGTCKPDKPNCKKS

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-Mlul




[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_007843

ORF Size: 207 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_007843.3
RefSeq Size:	896 bp
RefSeq ORF:	210 bp
Locus ID:	13214
UniProt ID:	P56386
Cytogenetics:	8 10.35 cM
Gene Summary:	Has bactericidal activity. May act as a ligand for C-C chemokine receptor CCR6. Positively regulates the sperm motility and bactericidal activity in a CCR6-dependent manner. Binds to CCR6 and triggers Ca ²⁺ mobilization in the sperm which is important for its motility. [UniProtKB/Swiss-Prot Function]