

Product datasheet for **MG200085**

Gng11 (NM_025331) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gng11 (NM_025331) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Gng11
Synonyms: 0610037B21Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200085 representing NM_025331
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTGCCCTTCACATCGAGGATCTGCCGAAAAGGAAAACTGAAGATGGAGTTGAGCAACTTCGCA
AAGAAGTCAAGTTCAGAGACAACAGGTATCTAAATGTTCTGAGGAAATAAAGAACTACATTGAAGAACG
TTCTGGAGAGGATCCTCTGGTAAAGGGAATCCAGAAGACAAGAATCCCTCAAAGAAAAGGCAGCTGT
GTCATTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200085 representing NM_025331
Red=Cloning site Green=Tags(s)
MPALHIEDLPEKEKLMKMEVEQLRKEVKLQRQQVSKCSEEIKNYIEERSGEDPLVKGIPEDKNPFKEKGSC
VIS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_025331

ORF Size: 219 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_025331.2 , NP_079607.1
RefSeq Size:	849 bp
RefSeq ORF:	222 bp
Locus ID:	66066
UniProt ID:	P61953
Cytogenetics:	6 A1
Gene Summary:	Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction (By similarity).[UniProtKB/Swiss-Prot Function]