

Product datasheet for **MG200111**

Gprc5c (BC008228) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gprc5c (BC008228) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Gprc5c
Synonyms: 1110028I06Rik; 3200002M13Rik; Raig3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200111 representing BC008228
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCACCGGAGGCACCCGATGGCCACCCAGGTCAGAGAGGGCCATTTTCAGGCCTGCAAGCTGCTGAAAC
 CAGTGGCTGAGATGGTGTTCATGGTACCAGACCCTGGACTCGAGATCTCCAGAACATCATTCTTTGTTCTA
 CATTTGGGTACCAGACAGAGCCTCTTCTATACCCTTTCCCCAAATGTGCTAACTGGCCCTCAGCCAG
 TCCTGCTTAGTGCCAAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200111 representing BC008228
 Red=Cloning site Green=Tags(s)
 MHRRHPMATQVREGHFQACKLLKPVAEMVSWYQTLDSRSPHHSLFYIWPDRASSIPFPNVLTWPSAQ
 SCLVPR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC008228

ORF Size: 230 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC008228
RefSeq Size:	1728 bp
RefSeq ORF:	230 bp
Locus ID:	70355
Cytogenetics:	11 E2
Gene Summary:	This retinoic acid-inducible G-protein coupled receptor provide evidence for a possible interaction between retinoid and G-protein signaling pathways.[UniProtKB/Swiss-Prot Function]