

## Product datasheet for **MG200128**

### Romo1 (NM\_025946) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Romo1 (NM\_025946) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Romo1  
**Synonyms:** 2010100O12Rik; AI853864  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG200128 representing NM\_025946  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCGGTGGCCGTGGTCCCTACGGGCAGTCCCAGCCAGCTGCTTCGACCGCGTGAAGATGGGCTTCG  
TCATGGGTTGCGCAGTGGGTATGGCGGCCGGGGCGCTGTTCCGGCACCTTCTCCTGTCTCAGGATCGGAAT  
GCGGGTGGGAGCTAATGGCGGCATTGGGAAAACCATGATGCAGAGTGGCGGCACGTTTGGCACTTTC  
ATGGCCATTGGAATGGGCATACGATGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG200128 representing NM\_025946  
Red=Cloning site Green=Tags(s)  
MPVAVGPYQSQPSCFDRVKMGFVMGCAVGMAGALFGTF SCLRIGMRGRELMGGIGKTMMSGGTFGTF  
MAIGMGIRC

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



[View online »](#)



<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_025946.6</a> , <a href="#">NP_080222.1</a>
<b>RefSeq Size:</b>	521 bp
<b>RefSeq ORF:</b>	240 bp
<b>Locus ID:</b>	67067
<b>UniProt ID:</b>	<a href="#">P60603</a>
<b>Cytogenetics:</b>	2 H1
<b>Gene Summary:</b>	Has antibacterial activity against a variety of bacteria including S.aureus, P.aeruginosa and M.tuberculosis. Acts by inducing bacterial membrane breakage (By similarity). [UniProtKB/Swiss-Prot Function]