

Product datasheet for **MG200134**

Cdc42se1 (NM_172395) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cdc42se1 (NM_172395) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Cdc42se1
Synonyms: 1300002M12Rik; AW558204; Cdcse1; SCIP1; Spec1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200134 representing NM_172395
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGAGCGAGTTTTGGCACAACTGGGCTGCTGCGTGGTCGAGAAGCCCCAGCCGAAGAAAAAGAGAAGAC
GGATTGACCGCACTATGATTGGGGAGCCAATGAACCTTGTCCACCTGACTCACATCGGCTCTGGGGAGAT
GGGGGCTGGAGATGGACTTGCCATGACAGGTGCAGTTCAGGAGCAGATGAGATCCAAGGGAACACCGT
GACAGACCGTGGAGCAATTCTAGGGCCTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200134 representing NM_172395
Red=Cloning site Green=Tags(s)

MSEFWHKLGCCVVEKPQPKKKRRRIDRTMIGEP MNFVHLTHIGSGEMGAGDGLAMTGAVQE QMRSGKNHR
 DRPWSNSRAL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI


[View online »](#)


Cloning Scheme:



ACCN:	NM_172395
ORF Size:	240 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_172395.3
RefSeq Size:	2757 bp

RefSeq ORF:	243 bp
Locus ID:	57912
UniProt ID:	Q8BHL7
Cytogenetics:	3 F2.1
Gene Summary:	Probably involved in the organization of the actin cytoskeleton by acting downstream of CDC42, inducing actin filament assembly. Alters CDC42-induced cell shape changes. In activated T-cells, may play a role in CDC42-mediated F-actin accumulation at the immunological synapse. May play a role in early contractile events in phagocytosis in macrophages (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG200134