

Product datasheet for **MG200186**

Sf3b5 (NM_175102) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: Sf3b5
Synonyms: 10kDa; I110005L13Rik; AU043053; Sf3b10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG200186 representing NM_175102
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACGGACCGGTACACCATCCACAGCCAGCTGGAGCATCTGCAGTCCAAGTACATCGGCACGGGCCACG
CCGACACCACCAAGTGGGAATGGCTGGTGAACCAGCACCGCGACTCCTACTGCTCCTACATGGGCCACTT
CGACCTCCTCAACTACTTCGCCATTGCCGAGAACGAGAGCAAAGCCCGCTGCGGTTCAACCTGATGGAG
AAAATGCTGCAGCCAGCGGGCCGCCCGCGGACAAGCCCGAGGAGAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200186 representing NM_175102
Red=Cloning site Green=Tags(s)

MTDRYTIHSQLEHLQSKYIGTGHADTTKWEVLVNQHRDSYCSYMGHFDLLNYFAIAENESKARVRFNLME
KMLQPSGPPADKPEEN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NM_175102

ORF Size: 258 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_175102.4 , NP_780311.2
RefSeq Size:	750 bp
RefSeq ORF:	261 bp
Locus ID:	66125
UniProt ID:	Q923D4
Cytogenetics:	10 A2

Gene Summary: Involved in pre-mRNA splicing as a component of the splicing factor SF3B complex, a constituent of the spliceosome. SF3B complex is required for 'A' complex assembly formed by the stable binding of U2 snRNP to the branchpoint sequence (BPS) in pre-mRNA. Sequence independent binding of SF3A/SF3B complex upstream of the branch site is essential, it may anchor U2 snRNP to the pre-mRNA.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG200186