

## Product datasheet for **MG200209**

### Atp5j2 (NM\_020582) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Atp5j2 (NM\_020582) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Atp5j2  
**Synonyms:** 1110019H14Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG200209 representing NM\_020582  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGTCGCTCGTGCCGCTGAAGGAGAAGAAGCTCATGGAAGTCAAACCTGGAGAGCTGCCGAGCTGGA  
TAATGATGCGGGATTTACCCCCAGTGGCATTGCAGGAGCCTTTCGGAGAGGGTATGACCGGTATTACAA  
CAAGTACATCAACGTTTCGAAAGGCAGCATCTCGGGATTAGCATGGTCTGGCAGCCTACGTGGTTTTTC  
AGCTACTGCATTTCTTACAAGGAACTCAAACATGAGCGGCGACGCAAGTACCAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG200209 representing NM\_020582  
Red=Cloning site Green=Tags(s)  
MASLVPLKEKKLMEVKLGELPSWIMMRDFTPSGIAGAFRRGYDRYYNKYINVRKGSISGISMVLAAYVVF  
SYCISYKELKHERRRKYH

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI




[View online »](#)

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_020582

**ORF Size:** 264 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_020582.1</a> , <a href="#">NP_065607.1</a>
<b>RefSeq Size:</b>	497 bp
<b>RefSeq ORF:</b>	267 bp
<b>Locus ID:</b>	57423
<b>UniProt ID:</b>	<a href="#">P56135</a>
<b>Cytogenetics:</b>	5 G2
<b>Gene Summary:</b>	Mitochondrial membrane ATP synthase (F(1)F(0) ATP synthase or Complex V) produces ATP from ADP in the presence of a proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain. F-type ATPases consist of two structural domains, F(1) - containing the extramembraneous catalytic core and F(0) - containing the membrane proton channel, linked together by a central stalk and a peripheral stalk. During catalysis, ATP synthesis in the catalytic domain of F(1) is coupled via a rotary mechanism of the central stalk subunits to proton translocation. Part of the complex F(0) domain. Minor subunit located with subunit a in the membrane.[UniProtKB/Swiss-Prot Function]