

Product datasheet for **MG200240**

Dpm3 (NM_026767) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dpm3 (NM_026767) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Dpm3
Synonyms: 1110001H19Rik; AV117483
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200240 representing NM_026767
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACGAAGTTAACACAGTGGCTTTGGGGACTGGCTCTCCTGGGCTCTGCCTGGGCTGCCCTGACCATGG
 GAGCACTGGGCTTGAGTTGCCTTTCCCTGCCGAGAGGTCCTGTGGCCACTGCCTGCCTACCTGTTGGT
 CTCCGCTGGCTGCTATGCCCTGGGCACGGTGGCTATCGCGTAGCTACATTTACGACTGCGAGGACGCT
 GCCCGAGAGCTGCAGAGCCAGATCGTGGAGCCCGAGCTGATTTAGCACGCAGGGGGCTGCGCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200240 representing NM_026767
 Red=Cloning site Green=Tags(s)
 MTKLTQWLWGLALLGSAWAALTMGALGLELPFPCREVLWPLPAYLLVSAGCYALGTVGYRVATFHDCEDA
 ARELQSQIVEARADLARRGLRF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_026767

ORF Size: 503 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026767.4 , NP_081043.1
RefSeq Size:	404 bp
RefSeq ORF:	279 bp
Locus ID:	68563
UniProt ID:	Q9D1Q4
Cytogenetics:	3 F1
Gene Summary:	Stabilizer subunit of the dolichol-phosphate mannose (DPM) synthase complex; tethers catalytic subunit DPM1 to the ER.[UniProtKB/Swiss-Prot Function]