

## Product datasheet for **MG200244**

### Ccl22 (NM\_009137) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccl22 (NM\_009137) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ccl22  
**Synonyms:** ABCD-1; DCBCK; MDC; Scya22  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG200244 representing NM\_009137  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTACCTGCGTGTCCCACTCCTGGTGGCTCTCGTCCTTCTTGCTGTGGCAATTCAGACCTCTGATG  
 CAGGTCCTATGGTGCCAATGTGGAAGACAGTATCTGCTGCCAGGACTACATCCGTCACCCTCTGCCATC  
 ACGTTTAGTGAAGGAGTTCTTCTGGACCTCAAAATCCTGCCGCAAGCCTGGCGTTGTTTTGATAACCGTC  
 AAGAACCAGATATCTGTGCCGATCCCAGGCAGGCTCTGGTGAAGAAGCTACTCCATAAACTGTCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG200244 representing NM\_009137  
 Red=Cloning site Green=Tags(s)  
 MATLRVPLLVALVLLAVAIQTS DAGPYGANVEDSICCDYIRHPLPSRLVKEFFWTSKSCRKPGVLLITV  
 KNRDICADPRQVWVKLLHKLS

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_009137

ORF Size: 276 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_009137.2</u>
<b>RefSeq Size:</b>	2225 bp
<b>RefSeq ORF:</b>	279 bp
<b>Locus ID:</b>	20299
<b>UniProt ID:</b>	<u>O88430</u>
<b>Cytogenetics:</b>	8 C5
<b>Gene Summary:</b>	Chemotactic for activated T-lymphocytes. May play an important role in the collaboration of dendritic cells and B-lymphocytes with T-cells in immune responses.[UniProtKB/Swiss-Prot Function]