

Product datasheet for **MG200303**

Prr15l (NM_146026) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prr15l (NM_146026) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Prr15l
Synonyms: Atad4; BC022765; E230006A12Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200303 representing NM_146026
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTGAAGTTGGTTGGTGAAGCTGACCTTCCTCCGAAAAAAAAATCCACTCCCAAAGTGCTCTACG
AAATCCCGGACACCTACGCCAAACAGAGGGAGGTGCAGAGCCCCAGGACCTGACGCAGGAGACCTCA
CAGCGACTTTAACAGCCGTCTGAAAAGATTGTGGATAAGAACACTAAAGGCAAGCATGTCAAGGTCTCC
AACTCAGGCCGTTCAAAGAGAAGAAAAAGTTCGGCCATGTTGGCAGAGAACCCCAACCTCTTTGATG
ACAGAGAGAACAAGGACAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200303 representing NM_146026
Red=Cloning site Green=Tags(s)
MTEVGWWKLTFLRKKKSTPKVLYEIPDTYAQTEGGAEPGPDAGDPHSDFNSRLEKIVDKNTKGKHKVKS
NSGRFKEKKKVRAMLAEPNLFDRENKGQ

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_146026

ORF Size: 300 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_146026.2
RefSeq Size:	1410 bp
RefSeq ORF:	303 bp
Locus ID:	217138
UniProt ID:	Q8JZM2
Cytogenetics:	11 D