

## Product datasheet for **MG200315**

### Ptms (NM\_026988) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Tag:** TurboGFP  
**Symbol:** Ptms  
**Synonyms:** 2610009E16Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)

**ORF Nucleotide Sequence:** >MG200315 representing NM\_026988  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTCGGAGAAGAGCGTGGAGGCAGCGCCGAGCTGAGTGCCAAGGACCTGAAAGAAAAGAAGGACAAGG  
TGGAGGAGAAGGCCGGCCGAAAGAACGGAAGAAAAGTAGTGGAGGAGGAGGAGAATGGAGCTGAGGA  
AGAGGAAGAAGAACTGCAGAGGATGGAGAGGATGATGATGAAGGGGATGAAGAAGATGAGGAAGAAGAG  
GAGGAGGATGAAGGCCCGTGCAGGAGAACTGCTGAAGAGGAGGATGAAGCAGACCCCAAGAGGCAGA  
AGACAGAAAACGGGCGTCGGCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG200315 representing NM\_026988  
Red=Cloning site Green=Tags(s)

MSEKSVEAAAELSAKDLKEKKDKVEEKAGRKERKKEVVEEEEENGAEIIIIIIETAEDEGDDDEGDEEIIIIII  
EEDEGPVVRKRTAEEDEADPKRQKTENGASA

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI



Cloning Scheme:



ACCN: NM\_026988

ORF Size: 303 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NM_026988.2</a>
<b>RefSeq Size:</b>	1184 bp
<b>RefSeq ORF:</b>	306 bp
<b>Locus ID:</b>	69202
<b>UniProt ID:</b>	<a href="#">Q9D0J8</a>
<b>Cytogenetics:</b>	6 F2
<b>Gene Summary:</b>	Parathyromosin may mediate immune function by blocking the effect of prothymosin alpha which confers resistance to certain opportunistic infections.[UniProtKB/Swiss-Prot Function]

**Product images:**

Circular map for MG200315