

Product datasheet for **MG200330**

Lsm3 (NM_026309) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lsm3 (NM_026309) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Lsm3
Synonyms: 1010001J12Rik; 2610005D18Rik; 6030401D18Rik; SMX4; USS2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200330 representing NM_026309
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGACGACGTAGATCAGCAACAGACCACCAATACCGTAGAGGAGCCCTGGATCTGATCAGGCTCA
GCCTGGATGAGCGAATTTATGTGAAAATGAGAAATGACCGAGAGCTTCGAGGCAGACTACACGCTTATGA
TCAACATTTAAATATGATCCTGGGAGATGTAGAAGAACTGTGACAACGATAGAGATTGATGAGGAGACA
TATGAAGAGATATATAAATCCACAAAACGAAACATCCCTATGCTCTTCGTCGGGGAGATGGCGTCGTT
TAGTCGCCCTCCATTGAGAGTTGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200330 representing NM_026309
Red=Cloning site Green=Tags(s)
MADDVDQQQTNTNTVEEPLDLIRLSLDERIYVKMRNDRELGRHLHAYDQHLNMILGDVEETVTTIEIDEET
YEEIYKSTKRNIPLFVRGDGVVLVAPPLRVG

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_026309

ORF Size: 306 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_026309.2
RefSeq Size:	634 bp
RefSeq ORF:	309 bp
Locus ID:	67678
UniProt ID:	P62311
Cytogenetics:	6 D1
Gene Summary:	Plays role in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex that is involved in spliceosome assembly, and as component of the precatalytic spliceosome (spliceosome B complex). The heptameric LSM2-8 complex binds specifically to the 3'-terminal U-tract of U6 snRNA.[UniProtKB/Swiss-Prot Function]