

Product datasheet for **MG200334**

Mycbp (NM_019660) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mycbp (NM_019660) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Mycbp
Synonyms: 5730488M09Rik; Amy-1; AW552132; AW742590
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200334 representing NM_019660
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCATTACAAAGCCCGCGACTCGAAGCGCGAGCAGTTCGGAGGTACTTGGAGAAGTCGGGGTGC
 TGGACACGCTGACGAAAGTGTGGTAGCCTTATGAAGAACCAGAGAAGCCACCAGTCTCTGGATTT
 TTTAAAGCATCACTTAGGAGCTGCTACCCCGAAAACCCAGAAAATAGAGCTGCTTCGCTAGAATTGGCA
 GAAATGAAAGAGAAATATGAAGCTACTGTAGAAGAAAATAAAAACTGAAAGCAAAGCTTGTTCAAGTATG
 AACCACTCAGGAGGAGAAGCGTCTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG200334 representing NM_019660
 Red=Cloning site Green=Tags(s)
 MAHYKAADSKREQFRYLEKSGVLDLTKVLVALYEEPEKPTSALDFLKHHLGAATPENPEIELLRLELA
 EMKEKYEATVEENKLLKAKLVQYEPPEEKRAE

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_019660

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019660.3</u>
RefSeq Size:	1486 bp
RefSeq ORF:	312 bp
Locus ID:	56309
UniProt ID:	<u>Q9EQS3</u>
Cytogenetics:	4 D2.2
Gene Summary:	May control the transcriptional activity of MYC. Stimulates the activation of E box-dependent transcription by MYC (By similarity).[UniProtKB/Swiss-Prot Function]