

Product datasheet for **MG200453**

Ltc4s (BC028760) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Ltc4s (BC028760) Mouse Tagged ORF Clone
 Tag: TurboGFP
 Symbol: Ltc4s
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >MG200453 representing BC028760
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGGACGAAGTGGCTCTTCTGGCTACCGTCACCCCTCGTGGGAGTTCTGTTGCAAGCCTACTTCTCCC
 TACAGGTGATCTCTGCACGAAGGGCTTCCACGTGTCGCCGCGCTCACCTCTGGCCCTCCCGAGTTCGA
 GCGCGTCTCCGAGCCAGGTAAGTGCAGCGAGTACTTCCGCTGTTCCCTCGCCACACTCTGGGTCGCC
 GGCATCTTCTCCACGAAGGAGCCGAGCCCTGTGCGGACTGTTCTACCTGTTGCGCGCCTCCGCTATT
 TCCAGGGATACGCGCGCTCAGCGCAACTCAGATGGCTCCAGATGCTCTGCCGATGGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200453 representing BC028760
 Red=Cloning site Green=Tags(s)

MKDEVALLATVTLVGVLLQAYFSLQVISARRAFHVSPPLTSGPPEFERVFRQVNCSEYFPLFLATLWVA
 GIFFHEGAAALCGLFYLFARLRYFGYARSAQLRWLQMLLPMA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC028760 , AAH28760
RefSeq Size:	503 bp
RefSeq ORF:	341 bp
Locus ID:	17001
Cytogenetics:	11 B1.3
Gene Summary:	The protein encoded by this gene is an enzyme that catalyzes the synthesis of leukotriene C4 by combining leukotriene A4 with reduced glutathione. The encoded protein is found in the outer nuclear membrane and in the peripheral endoplasmic reticulum. Leukotrienes have been implicated as mediators of anaphylaxis and inflammatory conditions such as bronchial asthma in humans. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]