

Product datasheet for **MG200460**

Rplp1 (BC092536) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rplp1 (BC092536) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rplp1
Synonyms: 2410042H16Rik; Arpp1; C430017H15Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200460 representing BC092536
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTCCGTCTCCGAGCTCGCTTGCATCTACTCCGCCCTCATCCTGCACGACGACGAGGTGACGGTCA
 CGGAGGATAAGATCAATGCCCTCATTAAAGCAGCTGGTGTGACGCGTGAACCTTTCTGGCCTGGCTTGT
 TGCCAAGGCTCTGGCCAATGTCAACATTGGGAGCCTCATCTGCAATGTAGGGGCTGGTGGCCCCGCTCCA
 GCAGCTGGAGCTGCGCTGCTGGTGGTGTCTCCATCCACTGCCGCCGCCAGCTGAGGAGAAGAAAG
 TGGAGGCAAAGAAGGAAGAGTCCGAGGAGTCTGAAGATGACATGGGCTTCGGTCTTTTGGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200460 representing BC092536
 Red=Cloning site Green=Tags(s)

MASVSELACIYSALILHDDEVTVTEKINALIKAAGVSVEPFWPGLFAKALANVNIGSLICNVGAGGPAP
 AAGAAPAGGAAPSTAAAPAEKKVEAKKEESESEDDMGFLFD

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC092536

ORF Size: 342 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC092536.1</u>
RefSeq Size:	627 bp
RefSeq ORF:	344 bp
Locus ID:	56040
Cytogenetics:	9 B
Gene Summary:	Plays an important role in the elongation step of protein synthesis.[UniProtKB/Swiss-Prot Function]