

## Product datasheet for **MG200554**

### Spata19 (BC049742) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Spata19 (BC049742) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Spata19  
**Synonyms:** 1700001L23Rik; AV047838; spergen-1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG200554 representing BC049742  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCGAGAGGATGTGGGGTTTTGGGCTGACGGAAGAGGAGGCTTCTCGGAGCATAAGGGAGAAGATGT  
CCATCAACGACTCTCCACCCATGGACATGACATACATGTGACCAGAGATTTGGTAAACACCATCTCTC  
AAAGTCTGATATGTTAACAGACCCGAGTCAAGAAGTCTGGAGGAAGAACAAGGATCCAGTTTATAAGA  
TGGAGCCACACCGTATCTTCCAAGTGCCAAGTGAAGTATGGATGACGTCATGCAGGAACGAATAGATC  
AAGTGAGACGAAGCGTATCTCATCTCATGTGTGACTCGTACAATGATCCAAGCTTCAGACTTCTTGCTC  
AGAATGC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG200554 representing BC049742  
Red=Cloning site Green=Tags(s)  
MSERMWGFGLTEEEASRSIREKMSINDSPHGHDIHVTRDLVKHHLKSDMLTDPSEQVLEERTRIQFIR  
WSHTRIFQVPSEVMDVMQERIDQVRRSVSHLMCDYNDPSFRTSCSEC

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC049742.1</a></u>
<b>RefSeq Size:</b>	671 bp
<b>RefSeq ORF:</b>	359 bp
<b>Locus ID:</b>	75469
<b>Cytogenetics:</b>	9 A4
<b>Gene Summary:</b>	May have a role in spermiogenesis.[UniProtKB/Swiss-Prot Function]