

Product datasheet for **MG200579**

Tmbim4 (BC027637) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmbim4 (BC027637) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Tmbim4
Synonyms: 0610007H07Rik; AU022431
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200579 representing BC027637
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGCCGACACGGACCCGGGCTATCCCCGCTCGTCCATCGAGGATGACTTCAACTACGGCAGCTGCGTGG
 CGTCGGCCAGCGTGACATCCGCATGGCCTTCTCAGAAAAGTCTACAGTATCCTCTCTGCAAGTCCT
 CCTGACTACAGTGACCTCTGCCCTGTTCTGTATTTCCAAGCTCTGCGGACATTTGTCCATGAAAGCCCT
 GCCTTAATTGTGGTGTCTGGGATCTCTGGGCTTGATCTTTGCACTGACTCTGCACAGACACACGC
 ATCCTCTGAACCTCTATCTACTCTTGCATTTGTAAGTACTTTAATACTGTCTACGCCATTAGCATATAA
 CCTGGGCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200579 representing BC027637
Red=Cloning site Green=Tags(s)

MADTDPGYPRSSIEDDFNYGSCVASASVHIRMAFLRKVYSILSLQVLLTTVTSLFLYFQALRTFVHESP
 ALIVVFALGSLGLIFALTLRHHTPLNLYLLFAFVSTLILSTPLAYNLGF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: BC027637

ORF Size: 360 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC027637](#), [AAH27637](#)

RefSeq Size: 504 bp

RefSeq ORF:	362 bp
Locus ID:	68212
Cytogenetics:	10 D2
Gene Summary:	Anti-apoptotic protein which can inhibit apoptosis induced by intrinsic and extrinsic apoptotic stimuli. Can modulate both capacitative Ca ²⁺ entry and inositol 1,4,5-trisphosphate (IP3)-mediated Ca ²⁺ release (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG200579