

Product datasheet for **MG200819**

Insl5 (NM_011831) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Insl5 (NM_011831) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Insl5
Synonyms: RIF2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200819 representing NM_011831
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGGGCCCCACTCTTGCTCTGTTTCTCTCTTAGTTCTGTTGGCTGTGGTGAAGTAAGAAGCAGGC
 AGACTGTGAAGCTCTGTGGCCTGGACTACGTGAGAACAGTTATCTACATCTGTCCAGCTCACGGTGGAG
 GAGACATCTGGAGGGCATTCCACTCTCAACAAGCTGAGACAAGAACTACCTCCAGCTCCTAGACAGG
 CACGAGCCATCCAAGAAACTCTGGAGCACAGCCTTCCAAGACGGATCTCTCAGGACAGGAGCTTGTT
 GAGATCCACAGGCACCCAAGGAAGTCTTTGGGAAGTGAAGAAGCACTCAGTGGTATCCAGACGAGATCT
 GCAAGCTCTGTGCTGCAGGGAAGGCTGCTCCATGAAGGAAGTCAAGCACTCTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200819 representing NM_011831
 Red=Cloning site Green=Tags(s)

MKGPTLALFLLLAVVEVRSRQTVKLCGLDYVRTVIYICASSRWRRLHLEGHFHSQQAETRNYLQLLDR
 HEPSKKTLEHSLPKTDLGQELVRDPQAPKEGLWELKKHSVVSRRDLQALCCREGCSMKELSTLC

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_011831

ORF Size: 405 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011831.2
RefSeq Size:	713 bp
RefSeq ORF:	408 bp
Locus ID:	23919
UniProt ID:	Q9WUG6
Cytogenetics:	4 47.27 cM
Gene Summary:	May have a role in gut contractility or in thymic development and regulation. Activates RXFP4 with high potency and appears to be the endogenous ligand for this receptor (By similarity). [UniProtKB/Swiss-Prot Function]