

Product datasheet for **MG200895**

Pfn1 (NM_011072) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pfn1 (NM_011072) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Pfn1
Synonyms: Pfn
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200895 representing NM_011072
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCGGGTGGAACGCCTACATCGACAGCCTTATGGCGGACGGGACCTGTCAGGACCGGCCATCGTAG
 GCTACAAGGACTCGCCCTCCGTCTGGGCCCGCGTCCCCGGGAAGACCTTCGTTAGCATTACGCCAGCTGA
 GGTTGGTGTCTGGTAGCAAAGACCGGTCAAGTTTTTCGTCAATGGGCTGACACTGGGGCCAGAAA
 TGTTCTGTGATCCGGGACTCACTGCTGCAAGACGGGAATTTACAATGGATCTTCGTACCAAGAGCACCG
 GAGGAGCCCCACCTTCAATGTCACTGTACCATGACTGCCAAGACGCTAGTCCTGCTGATGGGCAAAGA
 AGGTGTCCACGGTGGTTGATCAACAAGAAATGTTATGAAATGGCCTCTCACCTGCGGCGTCCCAGTAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200895 representing NM_011072
 Red=Cloning site Green=Tags(s)

MAGWNAYIDSLMADGTCQDAAIIVGYKDSPSVWAAVPGKTFVSIPTAEVGVLVGKDRSSFFVNGLTLGGQK
 CSVIRDSLLQDGEFTMDLRTKSTGGAPTFNVTVTMTAKTLVLLMGKEGVHGGLINKKCYEMASHLRRSQY

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_011072
 ORF Size: 420 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<p>NM_011072.4, NP_035202.1</p>
RefSeq Size:	<p>854 bp</p>
RefSeq ORF:	<p>423 bp</p>
Locus ID:	<p>18643</p>
UniProt ID:	<p>P62962</p>
Cytogenetics:	<p>11 43.21 cM</p>
Gene Summary:	<p>Binds to actin and affects the structure of the cytoskeleton. At high concentrations, profilin prevents the polymerization of actin, whereas it enhances it at low concentrations. By binding to PIP2, it inhibits the formation of IP3 and DG. Inhibits androgen receptor (AR) and HTT aggregation and binding of G-actin is essential for its inhibition of AR (By similarity). [UniProtKB/Swiss-Prot Function]</p>