

Product datasheet for **MG200954**

Eif1ax (NM_025437) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif1ax (NM_025437) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Eif1ax
Synonyms: 1500010B24Rik; AI426898; eIF-4C; Eif1ay
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG200954 representing NM_025437
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCCCAAGAATAAAGGAAAAGGAGGCAAAAACCGGCGTAGAGGCAAAATGAGAATGAATCTGAAAAGC
GAGAGCTTGTGTTCAAGGAGGATGGGCAAGAATATGCTCAAGTAATCAAATGTTGGAAACGGACGATT
AGAAGCCATGTGTTTGTGTTAAAGAGACTGTCCACATCAGAGGAAATTGAGAAAAAGGTTTG
ATAAATACCTCAGATATCATTTGGTGGGCTACGAGATTACCAGGACAACAAAGCCGATGTCATTTTAA
AATACAACGCAGATGAAGCTAGAAGTCTCAAGGCCTATGGCGAGCTTCCAGAGCATGCTAAAATCAATGA
AACTGATACATTTGGTCCTGGAGATGATGATGAAATTCAGTTTGATGACATTGGGGATGATGATGAAGAC
ATCGATGATATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG200954 representing NM_025437
Red=Cloning site Green=Tags(s)

MPKNKGKGGKNRRRGKNENESEKREL VKEDGQEYAQVIKMLGNRLEAMCFDGVKRLCHIRGKLRKKVW
 INTSDIILVGLRDYQDNKADVILKYNADERSLKAYGELPEHAKINETDTFGPGDDDEIQFDDIGDDDED
 IDDI


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_025437

ORF Size: 432 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025437.5</u>
RefSeq Size:	3019 bp
RefSeq ORF:	435 bp
Locus ID:	66235
UniProt ID:	<u>Q8BMJ3</u>
Cytogenetics:	X F4
Gene Summary:	Seems to be required for maximal rate of protein biosynthesis. Enhances ribosome dissociation into subunits and stabilizes the binding of the initiator Met-tRNA(I) to 40 S ribosomal subunits (By similarity).[UniProtKB/Swiss-Prot Function]