

Product datasheet for **MG201013**

Rpl27a (BC081430) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl27a (BC081430) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rpl27a
Synonyms: L27'; L27A; L29
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201013 representing BC081430
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGCCATCCAGACTGAGGAAGACCCGGAAACTCCGGGGCCACGTGAGCCACGGCCACGGCCGCATCGGTA
 AGCACCAGCAAGCAGCCAGGCGGCCGCGGAATGCTGGAGGCATGCACCACCACAGGATCAACTTTGACAA
 ATATCACCCAGGTTACTTTGGGAAAGTTGGTATGCGGCATTACCACTTGAAGAGGAACCAGAGCTTCTGC
 CCAACAGTCAACCTGGATAAACTGTGGACATTGGTCAGCGAGCAGACACGGGTCAATGCGGCAAAAAACA
 AGACTGGAGTTGCTCCCATCATTGATGTTGTTTCGATCAGGCTACTACAAAGTTCTGGCAAGGGAAAGCT
 CCCTAAGCAACCTGTATCGTGAAGCCAAATTCTTCAGCAGAAGAGCTGAGGAGAAGATTAAGGGTGT
 GGAGGTGCCTGTGTTCTGGTGGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201013 representing BC081430
 Red=Cloning site Green=Tags(s)

MPSRLRKTRKLRGHVSHGHGRIGKHKHPGGRNAGGMHHRINFDKYHPGYFGKVGMRHYHLKRNQSF
 PTVNLDKLVSEQTRVNAAKNKTGVAPIIDVVRSGYYKVLGKGLPKQPVIIVKAKFFSRRAEEKIKGV
 GGACVLVA

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC081430

ORF Size: 446 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC081430, AAH81430</u>
RefSeq Size:	518 bp
RefSeq ORF:	446 bp
Locus ID:	26451
Cytogenetics:	7 57.26 cM