

Product datasheet for **MG201028**

Calm4 (NM_020036) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Calm4 (NM_020036) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Calm4
Synonyms: 2310037J09Rik; A1649108; DD112; Scarf
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201028 representing NM_020036
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCACGGGTTTACTAAGGAGGAGGTCGCTGAGTTCAGGCAGCTTCAATAGGTTTGACAAGAATA
AGGATGGCCATATCAGCGTCGAGGAACCTGGAGATGTAATGAAGCAGCTGGCAAGAACCTCCAGAGAA
GGACCTGAAGGCTCTCATCTCCAAGTTGGACACAGATGGTGATGGCAAGATCAGCTTTGAAGAATCTTG
ACAGCCATAGAGAAGTACAAGAAAGGGCACAGAGCAGGGGAGCTGCGGGCTGTGTTCAATGTCCTTGACC
AGAATGGTGATGGCTACATCACTGTGGATGAGCTCAAGGAAAGCTTGCCAAGTTAGGCGAGTCCCTTTC
CCAGGAGGAGCTGGAGGACATGATCCGTGTGGCTGACGTGGACCAAGATGGGAAGGTGAAATATGAGGAG
TTTGTGCGCCTCCACATTGAAAAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201028 representing NM_020036
Red=Cloning site Green=Tags(s)

MSHGFTKEEVAEFQAAFNRFDKNKDGHISVEELGDVMQLGKNLPEKDLKALISKLDTDGDGKISFEEFL
TAIEKYKKGHRAGELRAVFNVLQNGDYITVDELKESLSKLGESLSQEELEDMIRVADVDQDGKVKYEE
FVRLHIEN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_020036

ORF Size: 444 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020036.2</u> , <u>NP_064420.2</u>
RefSeq Size:	916 bp
RefSeq ORF:	447 bp
Locus ID:	80796
UniProt ID:	<u>Q9JLM83</u>
Cytogenetics:	13 A1
Gene Summary:	Implicated in the early stage of ectopic ossification.[UniProtKB/Swiss-Prot Function]