

Product datasheet for **MG201096**

Rps18 (NM_011296) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rps18 (NM_011296) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rps18
Synonyms: H-2Ke3; H2-Ke3; Ke-3; ke3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201096 representing NM_011296
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTCTAGTGATCCCTGAGAAGTCCAGCACATTTTGGCAGTACTCAACACCAACATCGATGGCGGC
GGAAAATAGCCTTCGCCATCACTGCCATTAAGGGCGTGGGGCGGAGATATGCTCATGTGGTGTGAGGAA
AGCAGACATCGACCTACCAAGAGGGCTGGAGAACTCACGGAGGATGAGGTGGAGCGAGTGATCACCATC
ATGCAGAACCCACGACAGTACAAGATCCCAGACTGGTTCCTGAACAGACAGAAGGATGTGAAGGATGGGA
AGTACAGCCAGGTTCTGGCCAACGGTCTAGACAACAAGCTGCGTGAGGACCTGGAGAGGCTGAAGAAAAT
TCGAGCCCATAGAGGGCTGCGCCACTTTTGGGGCCTTCGTGTCCGCGGTCAGCACACCAAGACCCTGGC
CGCAGGGGCCGAAGTGTGGGTGTATCCAAGAAGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201096 representing NM_011296
Red=Cloning site Green=Tags(s)

MSLVIPEKFQHILRVLNTNIDGRRKIAFAITAIGVGRRYAHVVLRKADIDLTKRAGELTEDEVERVITI
MQNPRQYKIPDWFLNRQKDVKDGKYSQVLANGLDNKLRDLERLKKIRAHRLRHFVGLRVRGQHTKTTG
RRGRTVGVSKKK

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_011296

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_011296.3</u> |
| RefSeq Size: | 525 bp |
| RefSeq ORF: | 459 bp |
| Locus ID: | 20084 |
| UniProt ID: | <u>P62270</u> |
| Cytogenetics: | 17 17.98 cM |
| Gene Summary: | Located at the top of the head of the 40S subunit, it contacts several helices of the 18S rRNA. [UniProtKB/Swiss-Prot Function] |