

Product datasheet for **MG201120**

Ormdl1 (BC012315) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ormdl1 (BC012315) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Ormdl1
Synonyms: C730042F17Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201120 representing BC012315
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATGTTGGAGTTGCCACAGTGAGGTAACCCCAATACTCGTGAATGAACAGCCGTGGTATGTGGC
TGACATATGCGTTGGGAGTTGGCTGCTTCATATTGTATTACTCAGCATCCCCTTCTGCAGTGTTCTGT
TGCTGGACCTTAACCAACATAATACATAATCTGGGGATGTATGTGTTTTGCATGCAGTAAAGGAACG
CCTTTTGAGACTCCGGATCAGGGTAGAGCAAGGCTCCTAACTCACTGGGAACAGCTGGACTACGGAGTAC
AGTTTACATCTTCTCGGAAGTTTTTACAATTTCTCCAATAATTCTATACTTTCTGGCAAGTTTCTATAC
AAAGTATGACCCGACTCACTTCATCCTAAACACAGCTTCTCTGCTGAGTGTACTTATCCAAAATGCCA
CAGCTGCATGGTTCGCATCTTTGGAATTAATAAGTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201120 representing BC012315
Red=Cloning site **Green**=Tags(s)

MNVGVAHSEVNPNTVMNSRGMWLTALGVLLHIVLLSIPFCSVPVAVTLTNIHNLGMYVFLHAVKGT
PFETPDQGRARLLTHWEQLDYGVQFTSSRKFFITSPIILYFLASFYTKYDPTHFILNTASLLSVLIPKMP
QLHGVRIFGINKY

TRTRPLE - GFP Tag - V


Restriction Sites: Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: BC012315

ORF Size: 461 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC012315, AAH12315</u>
RefSeq Size:	1815 bp
RefSeq ORF:	461 bp
Locus ID:	227102
Cytogenetics:	1 C1.1
Gene Summary:	Negative regulator of sphingolipid synthesis.[UniProtKB/Swiss-Prot Function]