

Product datasheet for **MG201181**

Diablo (BC024780) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Diablo (BC024780) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Diablo
Synonyms: 0610041G12Rik; 1700006L01Rik; AU040403; Smac
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201181 representing BC024780
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGGCTCTGAGAAGTTGGGTGACTCGGAGCGTCTGCTCTCTGTTTCAGGTACAGACAGCGTTTTCTG
 TCTTGGCTAACTAAGAAACGCTGTTTCTCAGAATTGATAAAACCATGGCACAAAAGTGTGTTGACTGG
 ATTTGGCATGACTGTGTGCGTTCTATTGCTCAGAAATCGGAGCCTCAGTCTCTCAGTAACGAAGCA
 TTGATGAGGAGGGCTGTGCTTTGGTAACAGATAGCACCTCTACCTTTCTGTCTCAAACCACCTACGCGC
 TGATTGAAGCAATCACCGAGTATACTAAGGCTGTTTACACGTTAGTGTCTCTGTACCGACAATATACAAG
 TTTACTTGGGAAGATGAATCCAGGAGGAAGATGAGGTGTGGCAGGTAATTATAGGAGCCAGAGTTGAG
 GTGAGTGAACCTGTTAGTCATTTCCATATTAGTGGGCTTCATTATGGCATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201181 representing BC024780
 Red=Cloning site Green=Tags(s)

MAALRSWVTRSVCSLFRYRQRFVPLANSKKRCFSELIKPWKTVLTGFGMTLCAVPIAQKSEPQSL SNEA
 LMRRVSLVTDSTSTFLSQTTYALIEAITEYTKAVYTLVSLYRQYTSLLGKMNSQEEDVWQVIIGARVE
 VSEPVSHFHISGLHYGI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC024780

ORF Size: 471 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC024780.1</u>
RefSeq Size:	2585 bp
RefSeq ORF:	473 bp
Locus ID:	66593
Cytogenetics:	5 F
Gene Summary:	Promotes apoptosis by activating caspases in the cytochrome c/Apaf-1/caspase-9 pathway. Acts by opposing the inhibitory activity of inhibitor of apoptosis proteins (IAP). Inhibits the activity of BIRC6/bruce by inhibiting its binding to caspases (By similarity).[UniProtKB/Swiss-Prot Function]