

Product datasheet for **MG201192**

BC100379 Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: BC100379 Mouse Tagged ORF Clone
Tag: TurboGFP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG201192 representing BC100379
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTCTATGGGTCGGCAGGACCGCAGGCTGCAGCAGCAGCTTAATGAAAGCAGGAGCCGCTTCCAGA
CGTTAATGAAGCGGCTGATAGCGAAGTACAACCAGCCCTTAAAGGATAACTATGATGAGCTTAAAAAGA
ATTCAACAGACTTTATCAGAAGTATTGTCTGTACCGCAGCGGGCAAGGTGACTTCGTGTGCAGAGTA
TCGCCAATGAAAGCTGCTGCAGCTTCCCTGTCAATCAGAGCACTTGAAGAGACTAAATCCAGACTCCC
CACAAACAGAGTCCCAGAAACGGTCAATCTCTCCTGGCTGTACAGAAGAGTCTGCAGGATTCAACTCC
ACAGACTGCTAGTACACTGTGAGAGAGCCTTGGCTTCCAACAAAGAGCGCAAGTTGTCATACCGGTA
GCTTGTGCTCACCAGGCCAAATCTCACAATACCTCAGGAGCATCAGGCTGGCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG201192 representing BC100379
Red=Cloning site Green=Tags(s)

MESMGRQDRRLQQQLNESRSRFQTLMKRLIAKYNQPFKDNDELKKEFNRLYQKYCLSPQRAKVTSCGRV
SPMKAAAALPCQSEHLKRLNPDSPQSSQKRSISPGCHRRVLQDSTPQTASTLVREPWLPTKRRKLSYPV
ACAHQAKSHNTSGASGWP


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC100379

ORF Size: 474 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC100379.1
RefSeq Size:	2008 bp
RefSeq ORF:	476 bp