

## Product datasheet for **MG201285**

### 1600012H06Rik (BC027313) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 1600012H06Rik (BC027313) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** 1600012H06Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201285 representing BC027313  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTTTGAAGACGACGACGTCCTCCGGAGGAGTGGCTCCTTCTGCACGTGGTCCAGGGCCAAATAGGAGCCG  
 GGAATTACAGCTACTTGCAGTTGAACCATGAGGGTAAAATTATTCTCAGGATGCAGAGCCTGAGAGGGGA  
 TGCGGACCTGTACGTGTCGGACAGCACTCCGCACCCGAGCTTCGACGACTATGAGCTGCAGTCTGTCACT  
 TGCGGCCAGGATGTGGTCTCTATACCGGCACACTCCAGCGACCCGGTCCGGATAGGAATTTATGGACACC  
 CCTCTCACCATGAGAGCGACTTTGAGATGAGGGTATACTATGACAGAAGTGTGGATCAGTACCCATTCGG  
 AGAGGCTGCCTATTTCACTGACCCAACAGGAGCCAGTCAGCAGCAGGCTTACGCTCCAGAAGAGGCAGCC  
 CAAGAGGAGGAATCAGTTCTCTGGACGATATTGATTAGCATTGAAACTGGTGCTTGAATTCTCTTC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG201285 representing BC027313  
 Red=Cloning site Green=Tags(s)

MFEDDDVPEEWLLLHVQVQIGAGNYSYLRNLNHEGKIILRMQSLRGDADLYVSDSTPHPSFDDYELQSVT  
 CGQDVVSI PAHFQRPVGIYGHPSHHESDFEMRYYDRTVDQYPFGEAAFYTDPTGASQQAYAPEEAA  
 QEEESVLWTILISILKLVLEILF

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC027313

ORF Size: 491 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC027313</a></u> , <u><a href="#">AAH27313</a></u>
<b>RefSeq Size:</b>	2526 bp
<b>RefSeq ORF:</b>	491 bp
<b>Locus ID:</b>	67912
<b>Cytogenetics:</b>	17 A2
<b>Gene Summary:</b>	May be involved in induction of apoptosis in CD4(+) T-cells, but not CD8(+) T-cells or hepatocytes.[UniProtKB/Swiss-Prot Function]