

## Product datasheet for **MG201324**

### Rpl12 (BC081469) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rpl12 (BC081469) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Rpl12  
**Synonyms:** E430018F03  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG201324 representing BC081469  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCGCCAAGTTCGACCCCAACGAGGTCAAAGTCGTGTACCTGAGGTGCACCGGAGGCGAGGTCGGCG  
CCACATCCGCCTTGCCCCGAAGATCGGTCCTCTGGGTCTGTCTCCGAAGAAAGTTGGCGATGACATTGC  
CAAGGCTACCGGTGACTGGAAGGTCTCAGAATTACAGTGAAACTGACCATCCAGAACAGACAGGCCAG  
ATCGAGGTGGTGCCTTCTGCCTCAGCCCTGATCATCAAAGCCCTCAAGGAGCCACCGAGAGACAGGAAGA  
AGCAGAAGAACATTAAACACAGTGGGAACATTACTTTTGATGAGATTGTCAACATTGCCCGCAGATGAG  
ACACCGGTCTTTGGCTAGAGAACTTTCTGGAACATCAAGGAGATCCTGGGTACTGCACAGTCTGTGGGC  
TGCAATGTGGATGGCCGCCACCCTCATGACATCATAGATGACATCAACAGTGGTGCAGTGGAGTGCCAG  
CTAGT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG201324 representing BC081469  
**Red=Cloning site Green=Tags(s)**


MPPKFDPNVYVYLRCTGGEVGSALAPKIGPLGLSPKKVGGDDIAKATGDWKGLRITVKLTIQNRQAQ  
IEVVPASALIIKALKEPPDRKKQKNIKHSGNITFDEIVNIARQMRHRSRLARELSGTIKEILGTAQSVG  
CNVDGRHPHDIIDDINSGAVECPAS

**TRTRPLE** - GFP Tag - V


**Restriction Sites:** SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: BC081469

ORF Size: 497 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC081469</a></u> , <u><a href="#">AAH81469</a></u>
<b>RefSeq Size:</b>	640 bp
<b>RefSeq ORF:</b>	497 bp
<b>Locus ID:</b>	269261
<b>Cytogenetics:</b>	2 B
<b>Gene Summary:</b>	Binds directly to 26S ribosomal RNA.[UniProtKB/Swiss-Prot Function]